





OSFIX CATALOG PRODUCT

PT RISA IMPLANTAMA







ORTHOPAEDIC IMPLANTS

Orthopaedic implants are categorized into two groups including permanent joint replacements and temporary fracture fixation devices. Permanent orthopaedic implants include the shoulder, elbow, hip, knee, ankle, wrist, and finger joints, which are expected to serve in the human body throughout the life span of the patients. On the other hand, temporary orthopaedic implants including plates, screws, pins, wires, and intramedullary nails are needed to fix broken or fractured bones and are supposed to serve for a relatively short time just long enough to let bones heal. The types of implants that are installed vary depending on the location of the fracture and the fracture conditions that patient has suffered.

Properties of orthopaedic implants must ideally make the pieces chemically inert and absolutely biocompatible. Also, it must have a great strength, and a reasonable price so that it will be affordable by all. Another important feature is that it must be resist wearing for a long time.

BONE PLATES

OSFIX Implant Bone Plates are intended for repairing bone fractures, facilitating joint fusion, non – union and reconstruction of the orthopaedic extremities. The plates are anatomical shaped to fit the placement side. The bone plates are stainless steel plates which can accurately be adjusted to each body surface with the aid of the mentioned bending instruments. Bone Plates varies in thickness, length and form.

BONE SCREWS

OSFIX Implant Bone Screws includes different screw for repairing bone fracture, facilitating joint fusion, non – union and reconstruction of the orthopaedic extreimities. Bone screw are available as solid or cannulated.

BONE PIN

Bone pin such **OSFIX Schanz Screw and Steinmann Pin** are a metallic rod with pointed threads at one end and other end fashioned to fix into an attachment device. Schanz Screw self drilling has been specifically designed to optimise the bone / pin interface to reduce the occurrence of pin – tract related complications in external fixation. The unique design of Schanz Screw self drilling reduce heat generation and insertion torque while improving pullout resistance in cortical and cancellous bone.

ARTHROPLASTY (PROTHESIS)

Joint replacement is a procedure where by one or both the component forming a joint are replace with artificial component is called prosthesis. **OSFIX Prosthesis** are made up of special metal alloy to replace missing joints or bones, or to provide support to a damaged bone.



EXTERNAL FIXATION

Is a method of bony fixation outside the skin. With external fixation, a metal frame outside your body immobilizes the fracture with two or more pins that pass through your skin and into the bone on either side of the fracture. **OSFIX External Fixator** provides relative stability with little damage to soft tissue or disturbance to blood supply to bone. Simple, easy and minimally invasive stabilization.



1. MINI FRAGMENT PLATE

X Plate

Indication :

Used for mini fracture, fixation with 2.7 mm cortical screw



CODE	Holes
OF X 01	4
OF X 02	4

1.5 mm Mini T Plate A

Indication :

Used for fracture fixation of phalanges
Fixation with mini cortical screw



CODE	Holes
OF 10451.02	2
OF 10451.03	3
OF 10451.04	4
OF 10451.05	5

2.0 mm Mini T Plate B

Indication :

Used for fracture fixation of phalanges
Fixation with mini cortical screw



CODE	Holes
OF 10452.02	2
OF 10452.03	3
OF 10452.04	4
OF 10452.05	5



1.5 mm Mini T Plate C

Indication :

Used for fracture fixation of phalanges

Fixation with mini cortical screw



CODE	Holes
OF 10453.02	2
OF 10453.03	3
OF 10453.04	4
OF 10453.05	5

2.0 mm Mini T Plate D

Indication :

Used for fracture fixation of phalanges

Fixation with mini cortical screw



CODE	Holes
OF 10453.02	2
OF 10453.03	3
OF 10453.04	4
OF 10453.05	5

Mini T Oblique Plate (R/L)

Indication :

Used for fracture fixation of phalanges,
metatarsal, metacarpal

Fixation with mini cortical screw



CODE	Holes
OF 10471.02 R	2
OF 10471.02 L	2
OF 10471.03 R	3
OF 10471.03 L	3

OF 10471.04 R	4
OF 10471.04 L	4
OF 10471.05 R	5
OF 10471.05 L	5



1.5 mm Mini Straight Plate – DCP (Mini Compression Straight Plate)

Indication :

Used for fracture fixation of phalanges, metatarsal, zygoma, mid face skeleton and mandibular.



Fixation with mini cortical screw

Dynamic compression plate

CODE	Holes	Length (mm)
OF 1080.02	2	21
OF 1080.04	4	29
OF 1080.06	6	37

2.0 mm Mini Straight Plate – DCP (Mini Compression Straight Plate)

Indication :

Used for fracture fixation of phalanges, metatarsal, zygoma, mid face skeleton and mandibular.



Fixation with mini cortical screw

Dynamic compression plate

CODE	Holes	Length (mm)
OF 10801.02	2	21
OF 10801.04	4	29
OF 10801.06	6	37

1.5 mm Mini Curved Plate – DCP (Mini Compression Curved Plate)

Indication :

Used for fracture fixation of phalanges, metatarsal, zygoma, mid face skeleton and mandibular



Fixation with mini cortical screw

Dynamic compression plate



CODE	Holes
OF 1081.02	2
OF 1081.04	4
OF 1081.06	6

2.0 mm Mini Curved Plate – DCP (Mini Compression Curved Plate)

Indication :

Used for fracture fixation of phalanges, metatarsal, zygoma, mid face skeleton and mandibular

Fixation with mini cortical screw



Dynamic compression plate

CODE	Holes
OF 1081.02	2
OF 1081.04	4
OF 1081.06	6

Mini L Plate (R/L)

Indication :

Used for fracture fixation of phalanges, metatarsal, zygoma, mid face skeleton and mandibular

Fixation with mini cortical screw



CODE	Holes
OF 10480.02 R	2
OF 10480.02 L	2
OF 10480.03 R	3
OF 10480.03 L	3
OF 10480.04 R	4

OF 10480.04 L	4
OF 10480.05 R	5
OF 10480.05 L	5
OF 10480.06 R	6
OF 10480.06 L	6

Mini L Oblique Plate (R/L)

Indication :

Used for fracture fixation of phalanges, metatarsal, zygoma, mid face skeleton and mandibular

Fixation with mini cortical screw



CODE	Holes
OF 1149.02 R	2
OF 1149.02 L	2
OF 1149.03 R	3
OF 1149.03 L	3
OF 1149.04 R	4
OF 1149.04 L	4
OF 1149.05 R	5
OF 1149.05 L	5

1.5 mm Supermini Plate

Indication :

Used for cranio and severely comminuted fractures

Fixation with mini cortical screw



CODE	Holes
OF 10551.04	4
OF 10551.05	5
OF 10551.06	6
OF 10551.07	7
OF 10551.08	8
OF 10551.10	10

OF 10551.12	12
OF 10551.14	14
OF 10551.16	16
OF 10551.18	18
OF 10551.20	20
OF 10551.22	22
OF 10551.24	24

2.0 mm Supermini Plate

Indication :

Used for cranio and severely comminuted fractures

Fixation with mini cortical screw



CODE	Holes
OF 10552.04	4
OF 10552.05	5
OF 10552.06	6
OF 10552.07	7
OF 10552.08	8
OF 10552.10	10

OF 10552.12	12
OF 10552.14	14
OF 10552.16	16
OF 10552.18	18
OF 10552.20	20
OF 10552.22	22
OF 10552.24	24



Mini Condylar Plate

Indication :

Used for hand and foot mini condyle fractures

Fixation with mini cortical screw



CODE	Holes
OF 1031.03	3
OF 1031.04	4
OF 1031.05	5
OF 1031.06	6

Mini Compression T Plate

Indication :

Used for mid face skeleton and mandibular fractures

Fixation with mini cortical screw



CODE	Holes
OF 1090.02	2
OF 1090.04	4
OF 1090.06	6

Mini Compression L Plate (R/L)

Indication :

Used for mid face skeleton and mandibular fractures

Fixation with mini cortical screw



Available in Right and Left version

CODE	Holes
OF 1091.02	2
OF 1091.04	4
OF 1091.06	6



Mini Compression L Oblique Plate (R/L)

Indication :

Used for mid face skeleton and mandibular fractures

Fixation with mini cortical screw



Available in Right and Left version

CODE	Holes
OF 1092.02 R	2
OF 1092.02 L	2
OF 1092.04 R	4
OF 1092.04 L	4
OF 1092.06 R	6
OF 1092.06 L	6



2. MID FACE SKELETON PLATE

Maxillo Facial Plate

Indication :

Used for maxillo facial surgery



Fixation with mini cortical screw

CODE	Holes
OF 1040.04	4
OF 1040.05	5
OF 1040.06	6
OF 1040.07	7
OF 1040.08	8
OF 1040.09	9
OF 1040.10	10

OF 1040.11	11
OF 1040.12	12
OF 1040.14	14
OF 1040.16	16
OF 1040.18	18
OF 1040.20	20
OF 1040.22	22
OF 1040.24	24

Mini Maxillo Facial Plate

Indication :

Used for maxillo facial surgery



Fixation with mini cortical screw

CODE	Holes
OF 10401.04	4
OF 10401.05	5
OF 10401.06	6
OF 10401.07	7
OF 10401.08	8
OF 10401.09	9
OF 10401.10	10

OF 10401.11	11
OF 10401.12	12
OF 10401.14	14
OF 10401.16	16
OF 10401.18	18
OF 10401.20	20
OF 10401.22	22
OF 10401.24	24

Mandibula Plate

Indication :

Used for mandible surgery

Fixation with mini cortical screw



CODE	Holes
OF 4030.01	14
OF 4030.02	16
OF 4030.03	18

OF 4030.04	20
OF 4030.05	22
OF 4030.06	24



Angled Reconstruction Plate (R/L)

Indication :

Used for mandible surgery

Deep notches between holes

Fixation with 3.5 mm cortical screw



Available in Right and Left version

CODE	Holes
OF 800.23 R	23
OF 800.23 L	23
OF 800.24 R	24
OF 800.24 L	24

OF 800.28 R	28
OF 800.28 L	28
OF 800.31 R	31
OF 800.31 L	31

Angled Condyle Reconstruction Plate (R/L)

Indication :

Used for mandible surgery

Deep notches between holes

Fixation with 3.5 mm cortical screw



Available in Right and Left version

CODE	Holes
OF 810.23 R	23
OF 810.23 L	23
OF 810.24 R	24
OF 810.24 L	24
OF 810.28 R	28
OF 810.28 L	28
OF 810.31 R	31
OF 810.31 L	31



Bilateral Halfmoon Plate

Indication :

Used for mandible surgery

Deep notches between holes

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 810.23	23
OF 810.24	24
OF 810.28	28
OF 810.31	31

Mini Sherman Plate – Mandibula Reconstruction

Indication :

Used for mandible surgery

Fixation with mini cortical screw



CODE	Holes
OF 10701.04	4
OF 10701.06	6
OF 10701.08	8



3. SMALL FRAGMENT PLATE

1/3 Tubular Plate (One Third Tubular Plate)

Indication :

Used for fracture fixation of small bone, lateral malleolus, distal ulna, olecranon, distal humerus

Indicated for plating of fractures of long and small bone, the plate should only be used for load sharing purpose, e.g buttressing, tension bending, or neutralization

Oval holes – axial compression can be achieved

Fixation with 3.5 mm cortical screw or either 4.0 mm cancellous screw



CODE	Holes
OF 1030.03	3
OF 1030.04	4
OF 1030.05	5
OF 1030.06	6
OF 1030.07	7
OF 1030.08	8

OF 1030.09	9
OF 1030.10	10
OF 1030.11	11
OF 1030.12	12
OF 1030.13	13
OF 1030.14	14

1/4 Tubular Plate (Quarter Tubular Plate)

Indication :

Used for fracture fixation of small bone, metacarpal, metatarsal

Fixation with 2.7 mm cortical screw

Oval holes – axial compression can be achieved



CODE	Holes
OF 1031.03	3
OF 1031.04	4
OF 1031.05	5
OF 1031.06	6
OF 1031.07	7
OF 1031.08	8
OF 1031.09	9
OF 1031.10	10
OF 1031.11	11
OF 1031.12	12
OF 1031.13	13
OF 1031.14	14



3.5 mm Semi Tubular Plate

Indication :

Used for fracture fixation of small bone, lateral malleolus, distal ulna, olecranon, distal humerus

Indicated for plating of fractures of long and small bone, the plate should only be used for load sharing purpose, e.g buttressing, tension bending, or neutralization

Oval holes – axial compression can be achieved

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1010.04	4
OF 1010.05	5
OF 1010.06	6
OF 1010.07	7
OF 1010.08	8
OF 1010.09	9
OF 1010.10	10
OF 1010.11	11
OF 1010.12	12
OF 1010.13	13
OF 1010.14	14

Small Bone A Plate

Indication :

Used for fracture fixation of ulna and radius bone

Fixation with 3.5 mm cortical screw, the last two round holes on the fixation side can be used 4 mm cancellous screw

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



CODE	Holes
OF 1042.04	4
OF 1042.05	5
OF 1042.06	6
OF 1042.07	7
OF 1042.08	8

OF 1042.09	9
OF 1042.10	10
OF 1042.11	11
OF 1042.12	12
OF 1042.13	13
OF 1042.14	14



Small Bone A Plate - LC – AKD 21302810893

Indication :

Used for fracture fixation of ulna and radius bone

Fixation with 3.5 mm cortical screw, the last two round holes on the fixation side can be used 4 mm cancellous screw

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

CODE	Holes
OF 1042.04 LC	4
OF 1042.05 LC	5
OF 1042.06 LC	6
OF 1042.07 LC	7
OF 1042.08 LC	8

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



OF 1042.09 LC	9
OF 1042.10 LC	10
OF 1042.11 LC	11
OF 1042.12 LC	12
OF 1042.13 LC	13
OF 1042.14 LC	14

Small Bone A Plate - Locking

Indication :

Used for fracture fixation of ulna and radius bone

Indicated for plating of fractures of long and small bones, thr plate should only be used for load sharing purpose, e.g buttressing, tension banding, neutralization or compression

Fixation with 3.5 mm cortical screw, the last two round holes on the fixation side can be used 4 mm cancellous screw



CODE	Holes
OF 1042.04 L	4
OF 1042.05 L	5
OF 1042.06 L	6
OF 1042.07 L	7
OF 1042.08 L	8
OF 1042.09 L	9
OF 1042.10 L	10
OF 1042.11 L	11
OF 1042.12 L	12
OF 1042.13 L	13
OF 1042.14 L	14



Small Bone A Plate - LP

Indication :

Used for fracture fixation of ulna and radius bone

Fixation with 3.5 mm cortical screw, the last two round holes on the fixation side can be used 4 mm cancellous screw

Low profile plate facilitates fixation without impinging on soft tissue

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

Low profile plate

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



CODE	Holes
OF 1042.04 LP	4
OF 1042.05 LP	5
OF 1042.06 LP	6
OF 1042.07 LP	7
OF 1042.08 LP	8

OF 1042.09 LP	9
OF 1042.10 LP	10
OF 1042.11 LP	11
OF 1042.12 LP	12
OF 1042.13 LP	13
OF 1042.14 LP	14

Small Bone B Plate

Indication :

Used for fracture fixation of ulna and radius shaft (small adults and children)

Fixation with 3.5 mm cortical screw, the last two round holes on the fixation side can be used 4 mm cancellous screw



CODE	Holes
OF 1043.04	4
OF 1043.05	5
OF 1043.06	6
OF 1043.07	7
OF 1043.08	8
OF 1043.09	9
OF 1043.10	10
OF 1043.11	11
OF 1043.12	12
OF 1043.13	13
OF 1043.14	14



Small Bone C Plate

Indication :

Used for fracture fixation of ulna, radius, metatarsal and metacarpal bone

Fixation with 2.7 mm cortical screw



CODE	Holes
OF 1044.04	4
OF 1044.05	5
OF 1044.06	6
OF 1044.07	7
OF 1044.08	8
OF 1044.09	9
OF 1044.10	10
OF 1044.11	11
OF 1044.12	12
OF 1044.13	13
OF 1044.14	14

Small Bone D Plate

Indication :

Used for fracture fixation of ulna, and radius bone

Indicated for plating of fractures of long and small bone, the plate should only be used for load sharing purpose, e.g buttressing, tension bending, or neutralization

Fixation with 3.5 mm cortical screw or either 3.5 mm cancellous screw



Oval holes – axial compression can be achieved

CODE	Holes
OF 10421.04	4
OF 10421.05	5
OF 10421.06	6
OF 10421.07	7
OF 10421.08	8

OF 10421.09	9
OF 10421.10	10
OF 10421.11	11
OF 10421.12	12
OF 10421.13	13
OF 10421.14	14



T Oblique Plate (R/L)

Indication :

Used for fracture fixation of ulna, and radius bone

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw



Available in Right and Left version

CODE	Holes
OF 10461.02 R	2
OF 10461.02 L	2
OF 10461.03 R	3
OF 10461.03 L	3
OF 10461.04 R	4
OF 10461.04 L	4
OF 10421.05 R	5
OF 10421.05 L	5
OF 10421.06 R	6

OF 10421.06 L	6
OF 10421.07 R	7
OF 10421.07 L	7
OF 10421.08 R	8
OF 10421.08 L	8
OF 10421.09 R	9
OF 10421.09 L	9
OF 10421.10 R	10
OF 10421.10 L	10

T Oblique Plate - Locking (R/L)

Indication :

Used for fracture fixation of ulna, and radius bone

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw



Available in Right and Left version

CODE	Holes
OF 10461.02 L-R	2
OF 10461.02 L-L	2
OF 10461.03 L-R	3
OF 10461.03 L-L	3
OF 10461.04 L-R	4
OF 10461.04 L-L	4
OF 10421.05 L-R	5
OF 10421.05 L-L	5
OF 10421.06 L-R	6
OF 10421.06 L-L	6
OF 10421.07 L-R	7
OF 10421.07 L-L	7
OF 10421.08 L-R	8



OF 10421.08 L-L	8
OF 10421.09 L-R	9
OF 10421.09 L-L	9
OF 10421.10 L-R	10
OF 10421.10 L-L	10

Small T Plate Head 3

Indication :

Used for fracture fixation of distal radius bone

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw can be used in the head part of plate



CODE	Holes
OF 1046.02	2
OF 1046.03	3
OF 1046.04	4
OF 1046.05	5

OF 1046.06	6
OF 1046.07	7
OF 1046.08	8
OF 1046.09	9
OF. 1047.10	10

Small T Plate Head 3 – LP

Indication :

Used for fracture fixation of distal radius bone

Low profile plate facilitates fixation without impinging on soft tissue

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw can be used in the head part of plate

Low profile plate



CODE	Holes
OF 1046.02 LP	2
OF 1046.03 LP	3
OF 1046.04 LP	4
OF 1046.05 LP	5
OF 1046.06 LP	6
OF 1046.07 LP	7
OF 1046.08 LP	8
OF 1046.09 LP	9
OF. 1047.10 LP	10



Small T Plate Head 3 - LC

Indication :

Used for fracture fixation of distal radius, proximal and distal tibia, proximal humerus and clavicular

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw can be used in the head part of plate



Low contact plate

CODE	Holes
OF 1046.02 LC	2
OF 1046.03 LC	3
OF 1046.04 LC	4
OF 1046.05 LC	5
OF 1046.06 LC	6
OF 1046.07 LC	7
OF 1046.08 LC	8
OF 1046.09 LC	9
OF. 1047.10 LC	10

Small T Plate Head 3 – Locking

Indication :

Fracture fixation and fixation after osteotomies, malunions, non unions e.g including but not limited to distal radius, proximal and distal tibia, proximal humerus and clavicular

Fixation with 3.5 mm cortical locking screw and 4.0 mm cancellous locking screw can be used in the head part of plate



CODE	Holes
OF 1046.02 L	2
OF 1046.03 L	3
OF 1046.04 L	4
OF 1046.05 L	5
OF 1046.06 L	6
OF 1046.07 L	7
OF 1046.08 L	8
OF 1046.09 L	9
OF. 1047.10 L	10



Small T Plate Head 4

Indication :

Used for fracture fixation of distal radius, proximal and distal tibia, proximal humerus and clavicular

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw can be used in the head part of plate



CODE	Holes
OF 10464.02	2
OF 10464.03	3
OF 10464.04	4
OF 10464.05	5

OF 1046.06	6
OF 10464.07	7
OF 10464.08	8
OF 10464.09	9
OF. 10464.10	10

Small T Plate Head 4 – LP

Indication :

Used for fracture fixation of distal radius bone

Low profile plate facilitates fixation without impinging on soft tissue

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw can be used in the head part of plate

Low profile plate



CODE	Holes
OF 10464.02 LP	2
OF 10464.03 LP	3
OF 10464.04 LP	4
OF 10464.05 LP	5
OF 10464.06 LP	6
OF 10464.07 LP	7
OF 10464.08 LP	8
OF 10464.09 LP	9
OF. 10464.10 LP	10



Small T Plate Head 4 – LC

Indication :

Used for fracture fixation of distal radius bone

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw can be used in the head part of plate



Low contact plate

CODE	Holes
OF 10464.02 LC	2
OF 10464.03 LC	3
OF 10464.04 LC	4
OF 10464.05 LC	5
OF 10464.06 LC	6
OF 10464.07 LC	7
OF 10464.08 LC	8
OF 10464.09 LC	9
OF. 10464.10 LC	10

Small T Plate Head 4 – Locking

Indication :

Used for fracture fixation of distal radius bone

Fixation with 3.5 mm cortical locking screw and 4.0 mm cancellous locking screw can be used in the head part of plate



CODE	Holes
OF 10464.02 L	2
OF 10464.03 L	3
OF 10464.04 L	4
OF 10464.05 L	5
OF 10464.06 L	6
OF 10464.07 L	7
OF 10464.08 L	8
OF 10464.09 L	9
OF. 10464.10 L	10



Small T Buttress Plate

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 3.5 mm cortical screw and 4.0 mm cancellous screw



CODE	Holes
OF 10521.02 L	2
OF 10521.03 L	3
OF 10521.04 L	4

OF 10521.05 L	5
OF 10521.06 L	6
OF 10521.07 L	7
OF 10521.08 L	8
OF 10521.09 L	9
OF. 10521.10 L	10

Clavicula S Plate (R/L)

Indication :

Used for clavicle and calcaneus fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw

Available in Right and Left version



CODE	Holes
OF 1089.06 R	6
OF 1089.06 L	6
OF 1089.07 R	7
OF 1089.07 L	7
OF 1089.08 R	8
OF 1089.08 L	8
OF 1089.09 R	9

OF 1089.09 L	9
OF 1089.10R	10
OF 1089.10 L	10
OF 1089.11 R	11
OF 1089.11 L	11
OF 1089.12 R	12
OF 1089.12 L	12



Clavicula S Plate - Locking (R/L)

Indication :

Used for clavicle and calcaneus fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw



Available in Right and Left version

CODE	Holes
OF 1089.06 L-R	6
OF 1089.06 L-L	6
OF 1089.07 L-R	7
OF 1089.07 L-L	7
OF 1089.08 L-R	8
OF 1089.08 L-L	8
OF 1089.09 L-R	9

Locking plate

OF 1089.09 L-L	9
OF 1089.10 L-R	10
OF 1089.10 L-L	10
OF 1089.11 L-R	11
OF 1089.11 L-L	11
OF 1089.12 L-R	12
OF 1089.12 L-L	12

Clavicula Hook Plate (R/L)

Indication :

Fixation system with angular stability for lateral clavicle fractures and acromioclavicular joint injuries

Fixation with 3.5 mm cortical screw



Available in Right and Left version

CODE	Holes
OF 10892.03 R	3
OF 10892.03 L	3
OF 10892.04 R	4
OF 10892.04 L	4
OF 10892.05 R	5
OF 10892.05 L	5

OF 10892.06 R	6
OF 10892.06 L	6
OF 10892.07R	7
OF 10892.07 L	7
OF 10892.08 R	8
OF 10892.08 L	8



Clavicle Hook Plate - Locking (R/L)

Indication :

Fixation system with angular stability for lateral clavicle fractures and acromioclavicular joint injuries

Fixation with 3.5 mm cortical screw

Available in Right and Left version



CODE	Holes
OF 10892.03 L-R	3
OF 10892.03 L-L	3
OF 10892.04 L-R	4
OF 10892.04 L-L	4
OF 10892.05 L-R	5
OF 10892.05 L-L	5

OF 10892.06 L-R	6
OF 10892.06 L-L	6
OF 10892.07 L-R	7
OF 10892.07 L-L	7
OF 10892.08 L-R	8
OF 10892.08 L-L	8

Clavicle Z Plate (R/L)

Indication :

Fixation system with angular stability for lateral clavicle fractures and acromioclavicular joint injuries

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw

Available in Right and Left version



CODE	Holes
OF 10891.03 R	3
OF 10891.03 L	3
OF 10891.04 R	4
OF 10891.04 L	4
OF 10891.05 R	5
OF 10891.05 L	5

OF 10891.06 R	6
OF 10891.06 L	6
OF 10891.07R	7
OF 10891.07 L	7
OF 10891.08 R	8
OF 10891.08 L	8



Clover Leaf Plate

Indication :

With thinner head part, suitable for distal intra – articular tibial fractures and proximal humeral fractures

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw are used in the head part of plate



CODE	Holes
OF 1066.03 L	3
OF 1066.04 L	4
OF 1066.05 L	5
OF 1066.06 L	6
OF 1066.07 L	7

OF 1066.08 L	8
OF 1066.09 L	9
OF 1066.10 L	10
OF 1066.11 L	11
OF 1066.12 L	12

Clover Leaf Plate - Locking

Indication :

With thinner head part, suitable for distal intra – articular tibial fractures and proximal humeral fractures

For medial buttressing of comminuted distal tibial fractures

For comminuted humeral head fractures

Fixation with 3.5 mm cortical locking screw and 4 mm cancellous locking screw are used in the head part of plate



CODE	Holes
OF 1066.03 L	3
OF 1066.04 L	4
OF 1066.05 L	5
OF 1066.06 L	6
OF 1066.07 L	7
OF 1066.08 L	8
OF 1066.09 L	9
OF 1066.10 L	10
OF 1066.11 L	11
OF 1066.12 L	12



Clover Leaf Plate – LC

Indication :

With thinner head part, suitable for distal intra – articular tibial fractures and proximal humeral fractures

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw are used in the head part of plate



Low contact plate

CODE	Holes
OF 1066.03 LC	3
OF 1066.04 LC	4
OF 1066.05 LC	5
OF 1066.06 LC	6
OF 1066.07 LC	7

OF 1066.08 LC	8
OF 1066.09 LC	9
OF 1066.10 LC	10
OF 1066.11 LC	11
OF 1066.12 LC	12

Clover Leaf Plate – LP

Indication :

With thinner head part, suitable for distal intra – articular tibial fractures and proximal humeral fractures

Low profile plate facilitates fixation without impinging on soft tissue

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw are used in the head part of plate



Low profile plate

CODE	Holes
OF 1066.03 LP	3
OF 1066.04 LP	4
OF 1066.05 LP	5
OF 1066.06 LP	6
OF 1066.07 LP	7

OF 1066.08 LP	8
OF 1066.09 LP	9
OF 1066.10 LP	10
OF 1066.11 LP	11
OF 1066.12 LP	12



Calcaneal Plate (R/L)

Indication :

Used for calcaneus fractures

Fixation with 3.5 mm cortical screw and 3.5 mm cancellous screw



Available in Right and Left version

CODE	Holes
OF 1093.03 R	3
OF 1093.03 L	3
OF 1093.04 R	4
OF 1093.04 L	4
OF 1093.05 R	5
OF 1093.05 L	5

OF 1093.06 R	6
OF 1093.06 L	6
OF 1093.07R	7
OF 1093.07 L	7
OF 1093.08 R	8
OF 1093.08 L	8

Y Fork Plate (R/L)

Indication :

Used for distal humerus fractures

Fixation with 3.5 mm cortical screw and 3.5 cancellous screw

Available in Right and Left version



CODE	Holes
OF 10612.03 R	3
OF 10612.03 L	3
OF 10612.04 R	4
OF 10612.04 L	4
OF 10612.05 R	5
OF 10612.05 L	5
OF 10612.06 R	6
OF 10612.06 L	6
OF 10612.07 R	7
OF 10612.07 L	7

OF 10612.08 R	8
OF 10612.08 L	8
OF 10612.09 R	9
OF 10612.09 L	9
OF 10612.10 R	10
OF 10612.10 L	10
OF 10612.11 R	11
OF 10612.11 L	11
OF 10612.12 R	12
OF 10612.12 L	12



Condylar Humeral Plate – A (R/L)

Indication :

Used for distal humerus fractures

Fixation with 3.5 mm cortical screw and 3.5 mm cancellous screw



Available in Right and Left version

CODE	Holes
OF 1061.03 R	3
OF 1061.03 L	3
OF 1061.04 R	4
OF 1061.04 L	4
OF 1061.05 R	5
OF 1061.05 L	5
OF 1061.06 R	6
OF 1061.06 L	6
OF 1061.07 R	7
OF 1061.07 L	7

OF 1061.08 R	8
OF 1061.08 L	8
OF 1061.09 R	9
OF 1061.09 L	9
OF 1061.10 R	10
OF 1061.10 L	10
OF 1061.11 R	11
OF 1061.11 L	11
OF 1061.12 R	12
OF 1061.12 L	12

Condylar Humeral Plate – B (R/L)

Indication :

Used for distal humerus fractures

Fixation with 3.5 mm cortical screw and 3.5 mm cancellous screw



Available in Right and Left version

CODE	Holes
OF 10611.03 R	3
OF 10611.03 L	3
OF 10611.04 R	4
OF 10611.04 L	4
OF 10611.05 R	5
OF 10611.05 L	5
OF 10611.06 R	6
OF 10611.06 L	6
OF 10611.07 R	7
OF 10611.07 L	7

OF 10611.08 R	8
OF 10611.08 L	8
OF 10611.09 R	9
OF 10611.09 L	9
OF 10611.10 R	10
OF 10611.10 L	10
OF 10611.11 R	11
OF 10611.11 L	11
OF 10611.12 R	12
OF 10611.12 L	12



Humeral Plate (R/L)

Indication :

Used for proximal humerus, proximal tibia

Available in Right and Left version

Fixation with 3.5 mm cortical screws and 3.5 cancellous screws are used in the round holes



CODE	Holes
OF 10621.03 R	3
OF 10621.03 L	3
OF 10621.04 R	4
OF 10621.04 L	4
OF 10621.05 R	5
OF 10621.05 L	5
OF 10621.06 R	6
OF 10621.06 L	6

OF 10621.07 R	7
OF 10621.07 L	7
OF 10621.08 R	8
OF 10621.08 L	8
OF 10621.09 R	9
OF 10621.09 L	9
OF 10621.10 R	10
OF 10621.10 L	10

3.5 mm Reconstruction Plate – Straight

Indication :

Used for clavicle, calcaneus and pelvic acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw





CODE	Holes
OF 1068.04	4
OF 1068.05	5
OF 1068.06	6
OF 1068.07	7
OF 1068.08	8
OF 1068.09	9

OF 1068.10	10
OF 1068.11	11
OF 1068.12	12
OF 1068.13	13
OF 1068.14	14
OF 1068.15	15

3.5 mm Reconstruction Plate – Straight – Locking

Indication :

Used for clavicle, calcaneus and pelvic acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical locking screw



CODE	Holes
OF 1068.04 L	4
OF 1068.05 L	5
OF 1068.06 L	6
OF 1068.07 L	7
OF 1068.08 L	8
OF 1068.09 L	9

OF 1068.10 L	10
OF 1068.11 L	11
OF 1068.12 L	12
OF 1068.13 L	13
OF 1068.14 L	14
OF 1068.15 L	15

3.5 mm Reconstruction Plate – Curved

Indication :

Used for clavicle, calcaneus and pelvic acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1069.04	4
OF 1069.05	5
OF 1069.06	6
OF 1069.07	7
OF 1069.08	8
OF 1069.09	9

OF 1069.10	10
OF 1069.11	11
OF 1069.12	12
OF 1069.13	13
OF 1069.14	14
OF 1069.15	15



3.5 mm Reconstruction Plate – Curved – Locking

Indication :

Used for clavicle, calcaneus and pelvic acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical locking screw



CODE	Holes
OF 1069.04 L	4
OF 1069.05 L	5
OF 1069.06 L	6
OF 1069.07 L	7
OF 1069.08 L	8
OF 1069.09 L	9

OF 1069.10 L	10
OF 1069.11 L	11
OF 1069.12 L	12
OF 1069.13 L	13
OF 1069.14 L	14
OF 1069.15 L	15

Proximal Humeral Plate (R/L)

Indication :

Used for distal humerus fractures

Fixation with 3.5 mm cortical screw and 4 mm cancellous are used in the head part of plate

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments

Available in Right and Left version



CODE	Holes
OF 10622.03 R	3
OF 10622.03 L	3
OF 10622.04 R	4
OF 10622.04 L	4
OF 10622.05 R	5
OF 10622.05 L	5

OF 10622.06 R	6
OF 10622.06 L	6
OF 10622.07 R	7
OF 10622.07 L	7
OF 10622.08 R	8
OF 10622.08 L	8
OF 10622.09 R	9



OF 10622.09 L	9
OF 10622.10 R	10
OF 10622.10 L	10
OF 10622.11 R	11
OF 10622.11 L	11
OF 10622.12 R	12
OF 10622.12 L	12

OF 10622.13 R	13
OF 10622.13 L	13
OF 10622.14 R	14
OF 10622.14 L	14
OF 10622.15 R	15
OF 10622.15 L	15

Proximal Humeral Plate - Locking (R/L)

Indication :

Indicated for the osteosynthesis of fractures at tibia, humerus at proximal area

Fixation with 3.5 mm cortical locking screw and 4 mm cancellous screw are used in the head part of plate

Available in Right and Left version



CODE	Holes
OF 10622.03 L-R	3
OF 10622.03 L-L	3
OF 10622.04 L-R	4
OF 10622.04 L-L	4
OF 10622.05 L-R	5
OF 10622.05 L-L	5
OF 10622.06 L-R	6
OF 10622.06 L-L	6
OF 10622.07 L-R	7
OF 10622.07 L-L	7
OF 10622.08 L-R	8
OF 10622.08 L-L	8
OF 10622.09 L-R	9

OF 10622.09 L-L	9
OF 10622.10 L-R	10
OF 10622.10 L-L	10
OF 10622.11 L-R	11
OF 10622.11 L-L	11
OF 10622.12 L-R	12
OF 10622.12 L-L	12
OF 10622.13 L-R	13
OF 10622.13 L-L	13
OF 10622.14 L-R	14
OF 10622.14 L-L	14
OF 10622.15 L-R	15
OF 10622.15 L-L	15



Proximal Humeral Plate - LC (R/L)

Indication :

Used for distal humerus fractures

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw are used in the head part of plate

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments

Available in Right and Left version



CODE	Holes
OF 10622.03 R-LC	3
OF 10622.03 L-LC	3
OF 10622.04 R-LC	4
OF 10622.04 L-LC	4
OF 10622.05 R-LC	5
OF 10622.05 L-LC	5
OF 10622.06 R-LC	6
OF 10622.06 L-LC	6
OF 10622.07 R-LC	7
OF 10622.07 L-LC	7
OF 10622.08 R-LC	8
OF 10622.08 L-LC	8
OF 10622.09 R-LC	9

OF 10622.09 L-LC	9
OF 10622.10 R-LC	10
OF 10622.10 L-LC	10
OF 10622.11 R-LC	11
OF 10622.11 L-LC	11
OF 10622.12 R-LC	12
OF 10622.12 L-LC	12
OF 10622.13 R-LC	13
OF 10622.13 L-LC	13
OF 10622.14 R-LC	14
OF 10622.14 L-LC	14
OF 10622.15 R-LC	15
OF 10622.15 L-LC	15

Proximal Humeral Plate - LP (R/L)

Indication :

Used for distal humerus fractures

Low profile plate facilitates fixation without impinging on soft tissue

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw are used in the head part of plate

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder



The oval shape hole is used to create dynamic pressure between the fractured bone fragments

Available in Right and Left version

Low profile plate



CODE	Holes
OF 10622.03 R-LP	3
OF 10622.03 L-LP	3
OF 10622.04 R-LP	4
OF 10622.04 L-LP	4
OF 10622.05 R-LP	5
OF 10622.05 L-LP	5
OF 10622.06 R-LP	6
OF 10622.06 L-LP	6
OF 10622.07 R-LP	7
OF 10622.07 L-LP	7
OF 10622.08 R-LP	8
OF 10622.08 L-LP	8
OF 10622.09 R-LP	9

OF 10622.09 L-LP	9
OF 10622.10 R-LP	10
OF 10622.10 L-LP	10
OF 10622.11 R-LP	11
OF 10622.11 L-LP	11
OF 10622.12 R-LP	12
OF 10622.12 L-LP	12
OF 10622.13 R-LP	13
OF 10622.13 L-LP	13
OF 10622.14 R-LP	14
OF 10622.14 L-LP	14
OF 10622.15 R-LP	15
OF 10622.15 L-LP	15

Distal Medial Humeral Plate

Indication :

Used for distal medial humerus fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw and 3.5 mm cancellous screw





CODE	Holes
OF 1078.10	10
OF 1078.11	11
OF 1078.12	12
OF 1078.13	13
OF 1078.14	14
OF 1078.15	15
OF 1078.16	16
OF 1078.17	17
OF 1078.18	18
OF 1078.19	19
OF 1078.20	20

Distal Medial Humeral Plate - Locking

Indication :

Used for distal medial humerus fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical locking screw and 3.5 mm either 4mm cancellous locking screw

Locking plate



CODE	Holes
OF 1078.10 L	10
OF 1078.11 L	11
OF 1078.12 L	12
OF 1078.13 L	13
OF 1078.14 L	14
OF 1078.15 L	15
OF 1078.16 L	16
OF 1078.17 L	17
OF 1078.18 L	18
OF 1078.19 L	19
OF 1078.20 L	20



Fibula Plate

Indication :

Used for fibula fracture and small bone multiple (comminuted) fracture

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1063.07	7
OF 1063.08	8
OF 1063.09	9
OF 1063.10	10

OF 1063.11	11
OF 1063.12	12
OF 1063.13	13
OF 1063.14	14
OF 1063.15	15

3.5 mm Sherman Plate – Straight

Indication :

Used for pelvic acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1070.04	4
OF 1070.05	5
OF 1070.06	6
OF 1070.07	7
OF 1070.08	8
OF 1070.09	9
OF 1070.10	10
OF 1070.11	11
OF 1070.12	12
OF 1070.13	13
OF 1070.14	14
OF 1070.15	15

3.5 mm Sherman Plate – Straight – Locking

Indication :

Used for pelvis acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical locking screw





CODE	Holes
OF 1070.04 L	4
OF 1070.05 L	5
OF 1070.06 L	6
OF 1070.07 L	7
OF 1070.08 L	8
OF 1070.09 L	9

OF 1070.10 L	10
OF 1070.11 L	11
OF 1070.12 L	12
OF 1070.13 L	13
OF 1070.14 L	14
OF 1070.15 L	15

3.5 mm Sherman Plate – Curved

Indication :

Used for pelvis acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1071.04	4
OF 1071.05	5
OF 1071.06	6
OF 1071.07	7
OF 1071.08	8
OF 1071.09	9

OF 1071.10	10
OF 1071.11	11
OF 1071.12	12
OF 1071.13	13
OF 1071.14	14
OF 1071.15	15

3.5 mm Sherman Plate – Curved - Locking

Indication :

Used for pelvis acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 3.5 mm cortical locking screw

Locking plate



CODE	Holes
OF 1071.04 L	4
OF 1071.05 L	5
OF 1071.06 L	6
OF 1071.07 L	7
OF 1071.08 L	8
OF 1071.09 L	9

OF 1071.10 L	10
OF 1071.11 L	11
OF 1071.12 L	12
OF 1071.13 L	13
OF 1071.14 L	14
OF 1071.15 L	15



Slotted Narrow Plate

Indication :

Used for small bone fractures

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 200 87 A	2
OF 200 100 A	2
OF 200 107 A	2

OF 200 120 A	2
OF 200 150 A	2
OF 200 180 A	2
OF 200 200 A	2

3.5 mm Olecranon Hooked Plate

Indication :

Used for olecranon and proximal humerus fracture

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1067.04	4
OF 1067.05	5
OF 1067.06	6

OF 1067.07	7
OF 1067.08	8
OF 1067.09	9
OF 1067.10	10

Distal Humeral Y Plate (R/L)

Indication :

Used for intra articular fractures of distal humerus, supracondylar fractures, non unions of distal humerus

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw



Plate angled bending 40°

CODE	Holes
OF 1059.03 R	3
OF 1059.03 L	3
OF 1059.04 R	4
OF 1059.04 L	4
OF 1059.05 R	5
OF 1059.05 L	5

OF 1059.06 R	6
OF 1059.06 L	6
OF 1059.07 R	7
OF 1059.07 L	7
OF 1059.08 R	8
OF 1059.08 L	8



Rib Hooked Plate

Indication :

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw



CODE	Holes
OF 1083.12	3
OF 1083.14	4

OF 1083.16	5
OF 1083.18	6
OF 1083.20	7

Acromion Plate (R/L)

Indication :

Acromion fracture

Os acromiale

Scapular body fracture

Glenoid neck fracture

Intra – articular glenoid fracture

Scapula spine fracture

Fixation with 3.5 mm cortical screw and 4 mm cancellous screw



CODE	Holes
OF 1088.05 R	5
OF 1088.05 L	5
OF 1088.06R	6
OF 1088.06L	6
OF 1088.07R	7
OF 1088.07L	7




Distal Tibial Plate (R/L)

Indication :

Intended for fixation of complex intra and extra fractures and osteotomies of distal tibia

Fixation with 3.5 mm cortical screws, 4 mm cancellous screws are used in the head part of plate

Plate is twisted 20° and bent to fit distal tibia



CODE	Holes
OF 1073.04 R	4
OF 1073.04 L	4
OF 1073.05 R	5
OF 1073.05 L	5
OF 1073.06 R	6
OF 1073.06 L	6
OF 1073.07 R	7
OF 1073.07 L	7
OF 1073.08 R	8
OF 1073.08 L	8
OF 1073.09 R	9
OF 1073.09 L	9

OF 1073.10 R	10
OF 1073.10 L	10
OF 1073.11 R	11
OF 1073.11 L	11
OF 1073.12 R	12
OF 1073.12 L	12
OF 1073.13 R	13
OF 1073.13 L	13
OF 1073.14 R	14
OF 1073.14 L	14
OF 1073.15 R	15
OF 1073.15 L	15


Distal Tibial Plate (R/L) – LC

Indication :

Intended for fixation of complex intra and extra fractures and osteotomies of distal tibia

Fixation with 3.5 mm cortical screws, 4 mm cancellous screws are used in the head part of plate

Plate is twisted 20° and bent to fit distal tibia



CODE	Holes
OF 1073.04 R LC	4
OF 1073.04 L LC	4
OF 1073.05 R LC	5
OF 1073.05 L LC	5
OF 1073.06 R LC	6
OF 1073.06 L LC	6
OF 1073.07 R LC	7
OF 1073.07 L LC	7
OF 1073.08 R LC	8
OF 1073.08 L LC	8
OF 1073.09 R LC	9
OF 1073.09 L LC	9

OF 1073.10 R LC	10
OF 1073.10 L LC	10
OF 1073.11 R LC	11
OF 1073.11 L LC	11
OF 1073.12 R LC	12
OF 1073.12 L LC	12
OF 1073.13 R LC	13
OF 1073.13 L LC	13
OF 1073.14 R LC	14
OF 1073.14 L LC	14
OF 1073.15 R LC	15
OF 1073.15 L LC	15

Distal Tibial Plate (R/L) – Locking


Indication :

Intended for fixation of complex intra and extra fractures and osteotomies of distal tibia

The combination of dynamic compression hole with a locking screw hole provide the flexibility of axial compression and locking capability throughout the length of the plate shaft

Fixation with 3.5 mm cortical locking screws, 4 mm cancellous locking screws are used in the head part of plate

Plate is twisted 20° and bent to fit distal tibia



CODE	Holes
OF 1073.04 LR	4
OF 1073.04 LL	4
OF 1073.05 LR	5
OF 1073.05 LL	5
OF 1073.06 LR	6
OF 1073.06 LL	6
OF 1073.07 LR	7
OF 1073.07 LL	7
OF 1073.08 LR	8
OF 1073.08 LL	8
OF 1073.09 LR	9
OF 1073.09 LL	9

OF 1073.10 LR	10
OF 1073.10 LL	10
OF 1073.11 LR	11
OF 1073.11 LL	11
OF 1073.12 LR	12
OF 1073.12 LL	12
OF 1073.13 LR	13
OF 1073.13 LL	13
OF 1073.14 LR	14
OF 1073.14 LL	14
OF 1073.15 LR	15
OF 1073.15 LL	15



4. LARGE FRAGMENT PLATE

4.5 mm Semi Tubular Plate

Indication :

Used for fracture fixation of the forearm, subcutaneous of tibia, neutralization plate for radius and fibula

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw



CODE	Holes
OF 1010.04	4
OF 1010.05	5
OF 1010.06	6
OF 1010.07	7
OF 1010.08	8
OF 1010.09	9
OF 1010.10	10
OF 1010.11	11
OF 1010.12	12
OF 1010.13	13

Narrow Plate – ACP

Indication :

As neutralization and tension band plate on tibia and ulna (humerus, forearm)

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

The screw holes in the ACP are the best describes as a portion of an inclined and angled cylinder

Like a ball, the screw head slide down the inclined shoulder of the cylinder.

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



CODE	Holes
OF 1020.04	4
OF 1020.05	5
OF 1020.06	6
OF 1020.07	7
OF 1020.08	8

OF 1020.09	9
OF 1020.10	10
OF 1020.11	11
OF 1020.12	12
OF 1020.13	13
OF 1020.14	14



OF 1020.15	15
OF 1020.16	16
OF 1020.17	17

OF 1020.18	18
OF 1020.19	19
OF 1020.20	20

Narrow Plate ACP – LC

Indication :

As neutralization and tension band plate on tibia and ulna (humerus, forearm)

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

Like a ball, the screw head slide down the inclined shoulder of the cylinder.

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



CODE	Holes
OF 1020.04 LC	4
OF 1020.05 LC	5
OF 1020.06 LC	6
OF 1020.07 LC	7
OF 1020.08 LC	8
OF 1020.09 LC	9
OF 1020.10 LC	10
OF 1020.11 LC	11

OF 1020.12 LC	12
OF 1020.13 LC	13
OF 1020.14 LC	14
OF 1020.15 LC	15
OF 1020.16 LC	16
OF 1020.17 LC	17
OF 1020.18 LC	18
OF 1020.19 LC	19
OF 1020.20 LC	20

Narrow Plate ACP – LP

Indication :

As neutralization and tension band plate on tibia and ulna (humerus, forearm)

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Low profile plate facilitates fixation without impinging on soft tissue

The screw holes in the ACP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments





Low profile plate

CODE	Holes
OF 1020.04 LP	4
OF 1020.05 LP	5
OF 1020.06 LP	6
OF 1020.07 LP	7
OF 1020.08 LP	8
OF 1020.09 LP	9
OF 1020.10 LP	10
OF 1020.11 LP	11

OF 1020.12 LP	12
OF 1020.13 LP	13
OF 1020.14 LP	14
OF 1020.15 LP	15
OF 1020.16 LP	16
OF 1020.17 LP	17
OF 1020.18 LP	18
OF 1020.19 LP	19
OF 1020.20 LP	20

Narrow Plate ACP – Locking

Indication :

As neutralization and tension band plate on tibia and ulna (humerus, forearm)

For extra – articular fractures of the metaphyseal area that extend into the shaft area

Also use in fixation of periprosthetic fractures, osteopenic bone and non unions or malunions

Locking compression plate are intended for temporary fixation, correction or stabilization of bones in various anatomical regions

Fixation with 4.5 mm cortical locking screw or either 6.5 mm cancellous locking screw



CODE	Holes
OF 1020.04 L	4
OF 1020.05 L	5
OF 1020.06 L	6
OF 1020.07 L	7
OF 1020.08 L	8
OF 1020.09 L	9
OF 1020.10 L	10
OF 1020.11 L	11

OF 1020.12 L	12
OF 1020.13 L	13
OF 1020.14 L	14
OF 1020.15 L	15
OF 1020.16 L	16
OF 1020.17 L	17
OF 1020.18 L	18
OF 1020.19 L	19
OF 1020.20 L	20

Narrow Plate – Round Hole

Indication :

As neutralization and tension band plate on tibia and ulna (humerus, forearm)

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw





CODE	Holes
OF 1020.04R	4
OF 1020.05R	5
OF 1020.06R	6
OF 1020.07R	7
OF 1020.08R	8
OF 1020.09R	9
OF 1020.10R	10
OF 1020.11R	11

OF 1020.12R	12
OF 1020.13R	13
OF 1020.14R	14
OF 1020.15R	15
OF 1020.16R	16
OF 1020.17R	17
OF 1020.18R	18
OF 1020.19R	19
OF 1020.20R	20

Narrow Plate – Bridging

Indication :

Used for tumor resection, bridging or bone elongation

Instead of individually fixing each fragment

Minimal disruption to blood supply

Reduction is performed indirectly

Hole distance : 12 mm

CODE	Holes	Elongation (mm)
OF 1020.06 B	6	30
OF 1020.06 B	6	40
OF 1020.06 B	6	50
OF 1020.06 B	6	60
OF 1020.06 B	6	80
OF 1020.08 B	8	30
OF 1020.08 B	8	40
OF 1020.08 B	8	50
OF 1020.08 B	8	60
OF 1020.08 B	8	80
OF 1020.10 B	10	30
OF 1020.10 B	10	40
OF 1020.10 B	10	50
OF 1020.10 B	10	60
OF 1020.10 B	10	80



Broad / Wide Plate – ACP

Indication :

As tension band plate on femur and for humeral pseudarthrosis

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



Auto compression plate

CODE	Holes
OF 1021.04	4
OF 1021.05	5
OF 1021.06	6
OF 1021.07	7
OF 1021.08	8
OF 1021.09	9
OF 1021.10	10
OF 1021.11	11

OF 1021.12	12
OF 1021.13	13
OF 1021.14	14
OF 1021.15	15
OF 1021.16	16
OF 1021.17	17
OF 1021.18	18
OF 1021.19	19
OF 1021.20	20

Broad / Wide Plate - ACP – LC

Indication :

As tension band plate on femur and for humeral pseudarthrosis

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

The oval shape hole is used to create dynamic pressure between the fractured bone fragments





Low contact plate

CODE	Holes
OF 1021.04 LC	4
OF 1021.05 LC	5
OF 1021.06 LC	6
OF 1021.07 LC	7
OF 1021.08 LC	8
OF 1021.09 LC	9
OF 1021.10 LC	10
OF 1021.11 LC	11

OF 1021.12 LC	12
OF 1021.13 LC	13
OF 1021.14 LC	14
OF 1021.15 LC	15
OF 1021.16 LC	16
OF 1021.17 LC	17
OF 1021.18 LC	18
OF 1021.19 LC	19
OF 1021.20 LC	20

Broad / Wide Plate – ACP – LP

Indication :

As tension band plate on femur and for humeral pseudarthrosis

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Low profile plate facilitates fixation without impinging on soft tissue

The screw holes in the DCP are the best describes as a portion of an inclined and angled cylinder

Low profile plate

CODE	Holes
OF 1021.04 LP	4
OF 1021.05 LP	5
OF 1021.06 LP	6
OF 1021.07 LP	7
OF 1021.08 LP	8
OF 1021.09 LP	9
OF 1021.10 LP	10
OF 1021.11 LP	11

The oval shape hole is used to create dynamic pressure between the fractured bone fragments



OF 1021.12 LP	12
OF 1021.13 LP	13
OF 1021.14 LP	14
OF 1021.15 LP	15
OF 1021.16 LP	16
OF 1021.17 LP	17
OF 1021.18 LP	18
OF 1021.19 LP	19
OF 1021.20 LP	20

Broad / Wide Plate – Locking

Indication :

For extra – articular fractures of the metaphyseal area that extend into the shaft area

For the osteosynthesis of fractures at femur and tibia at shaft area

Also use in fixation of periprosthetic fractures, osteopenic bone and non unions or malunions



Fixation with 4.5 mm cortical locking screw or either 6.5 mm cancellous locking screw

Locking compression plate are intended for temporary fixation, correction or stabilization of bones in various anatomical regions



CODE	Holes
OF 1021.04 L	4
OF 1021.05 L	5
OF 1021.06 L	6
OF 1021.07 L	7
OF 1021.08 L	8
OF 1021.09 L	9
OF 1021.10 L	10
OF 1021.11 L	11
OF 1021.12 L	12
OF 1021.13 L	13

OF 1021.14 L	14
OF 1021.15 L	15
OF 1021.16 L	16
OF 1021.17 L	17
OF 1021.18 L	18
OF 1021.19 L	19
OF 1021.20 L	20

Broad / Wide Plate – Round Hole

Indication :

As tension band plate on femur and for humeral pseudarthrosis

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw



CODE	Holes
OF 1021.04 R	4
OF 1021.05 R	5
OF 1021.06 R	6
OF 1021.07 R	7
OF 1021.08 R	8
OF 1021.09 R	9
OF 1021.10 R	10
OF 1021.11 R	11
OF 1021.12 R	12

OF 1021.13 R	13
OF 1021.14 R	14
OF 1021.15 R	15
OF 1021.16 R	16
OF 1021.17 R	17
OF 1021.18 R	18
OF 1021.19 R	19
OF 1021.20 R	20



Broad / Wide Plate – Bridging

Indication :

Tumor resection, bridging or bone elongation

Instead of individually fixxing each fragment

Minimal disruption to blood supply

Reduction is performed indirectly

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw



Hole distance : 12 mm

Elongation : 30; 40; 50; 60 mm

CODE	Holes	Elongation (mm)
OF 1021.06 B	6	30
OF 1021.06 B	6	40
OF 1021.06 B	6	50
OF 1021.06 B	6	60
OF 1021.06 B	6	80
OF 1021.08 B	8	30
OF 1021.08 B	8	40

OF 1021.08 B	8	50
OF 1021.08 B	8	60
OF 1021.08 B	8	80
OF 1021.10 B	10	30
OF 1021.10 B	10	40
OF 1021.10 B	10	50
OF 1021.10 B	10	60
OF 1021.10 B	10	80

Slotted Broad Plate

Indication :

Fixation with 4.5 mm cortical screw

Used for large bone fractures



CODE	Holes (Slot)
OF 200 87 B	2
OF 200 100 B	2
OF 200 107 B	2

OF 200 120 B	2
OF 200 150 B	2
OF 200180 B	2
OF 200200 B	2

Egger Plate

Indication :

Used for large bone fractures

Fixation with 4.5 mm cortical screw



CODE	Holes (Slot)
OF 210 87	4
OF 210 100	4
OF 210 107	4
OF 210 120	4
OF 210 150	4
OF 210 180	4
OF 210 200	4

T Plate

Indication :

Shaft with round profile for humeral and tibial head

Intended to buttress metaphyseal fractures of proximal humerus, media tibial plateau and distal tibia.

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw



CODE	Holes
OF 1051.03	3
OF 1051.04	4
OF 1051.05	5
OF 1051.06	6
OF 1051.07	7
OF 1051.08	8

OF 1051.09	9
OF 1051.10	10
OF 1051.11	11
OF 1051.12	12
OF 1051.13	13
OF 1051.14	14
OF 1051.15	15

T Plate - Locking

Indication :

Intended to buttress metaphyseal fractures of proximal humerus, media tibial plateau and distal tibia.

Also for use in fixation of osteopenic bone and fixation of non unions and malunions

Fixation with 4.5 mm cortical locking screw or either 6.5 mm cancellous locking screw



CODE	Holes
OF 1051.03 L	3
OF 1051.04 L	4
OF 1051.05 L	5
OF 1051.06 L	6
OF 1051.07 L	7
OF 1051.08 L	8

OF 1051.09 L	9
OF 1051.10 L	10
OF 1051.11 L	11
OF 1051.12 L	12
OF 1051.13 L	13
OF 1051.14 L	14
OF 1051.15 L	15

T Plate - LP

Indication :

Shaft with raound profile for humeral and tibial head

Low profile plate facilitates fixation without impinging on soft tissue

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw



Low profile plate

CODE	Holes
OF 1051.03 LP	3
OF 1051.04 LP	4
OF 1051.05 LP	5
OF 1051.06 LP	6
OF 1051.07 LP	7

OF 1051.08 LP	8
OF 1051.09 LP	9
OF 1051.10 LP	10
OF 1051.11 LP	11
OF 1051.12 LP	12
OF 1051.13 LP	13
OF 1051.14 LP	14
OF 1051.15 LP	15



Form Plate (R/L)

Indication :

Used for distal tibia

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head of plate

Available in Right and Left version



CODE	Holes
OF 10512.03 R	3
OF 10512.03 L	3
OF 10512.04 R	4
OF 10512.04 L	4
OF 10512.05 R	5
OF 10512.05 L	5
OF 10512.06 R	6
OF 10512.06 L	6

OF 10512.07 R	7
OF 10512.07 L	7
OF 10512.08 R	8
OF 10512.08 L	8
OF 10512.09 R	9
OF 10512.09 L	9
OF 10512.10 R	10
OF 10512.10 L	10

Tibia L Plate (R/L)

Indication :

Used for tibia fractures

Fixation with 4.5 mm cortical screws or either 6.5 mm cancellous screws

Available in Right and Left version

Plate angled bending 20°



CODE	Holes
OF 10513.03 R	3
OF 10513.03 L	3
OF 10513.04 R	4
OF 10513.04 L	4
OF 10513.05 R	5
OF 10513.05 L	5
OF 10513.06 R	6
OF 10513.06 L	6

OF 10513.07 R	7
OF 10513.07 L	7
OF 10513.08 R	8
OF 10513.08 L	8
OF 10513.09 R	9
OF 10513.09 L	9
OF 10513.10 R	10
OF 10513.10 L	10



Tibia L Plate (R/L) – Locking

Indication :

Used for proximal tibial fractures

Fixation with 4.5 mm cortical locking screws or either 6.5 mm cancellous locking screws

Available in Right and Left version

Plate angled bending 20°



CODE	Holes
OF 10513.03 L-R	3
OF 10513.03 L-L	3
OF 10513.04 R-L	4
OF 10513.04 L-L	4
OF 10513.05 L-R	5
OF 10513.05 L-L	5
OF 10513.06 L-R	6
OF 10513.06 L-L	6

OF 10513.07 L-R	7
OF 10513.07 L-L	7
OF 10513.08 L-R	8
OF 10513.08 L-L	8
OF 10513.09 L-R	9
OF 10513.09 L-L	9
OF 10513.10 L-R	10
OF 10513.10 L-L	10

T Buttress Plate

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 4.5 mm cortical screws or either 6.5 mm cancellous screws



CODE	Holes
OF 1052.03	3
OF 1052.04	4
OF 1052.05	5
OF 1052.06	6

OF 1052.07	7
OF 1052.08	8
OF 1052.09	9
OF 1052.10	10



T Buttress Plate – Locking

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 4.5 mm cortical locking screws or either 6.5 mm cancellous locking screws



CODE	Holes
OF 1052.03 L	3
OF 1052.04 L	4
OF 1052.05 L	5
OF 1052.06 L	6
OF 1052.07 L	7
OF 1052.08 L	8
OF 1052.09 L	9
OF 1052.10 L	10

L Buttress Plate (R/L)

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 4.5 mm cortical screws or either 6.5 mm cancellous screws



Available in Right and Left version

CODE	Holes
OF 1053.03 R	3
OF 1053.03 L	3
OF 1053.04 R	4
OF 1053.04 L	4
OF 1053.05 R	5

OF 1053.05 L	5
OF 1053.06 R	6
OF 1053.06 L	6
OF 1053.07 R	7
OF 1053.07 L	7
OF 1053.08 R	8



OF 1053.08 L	8
OF 1053.09 R	9
OF 1053.09 L	9

OF 1053.10 R	10
OF 1053.10 L	10

L Buttress Plate (R/L) – Locking

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Available in Right and Left version

CODE	Holes
OF 1053.03 L-R	3
OF 1053.03 L-L	3
OF 1053.04 L-R	4
OF 1053.04 L-L	4
OF 1053.05 L-R	5
OF 1053.05 L-L	5
OF 1053.06 L-R	6
OF 1053.06 L-L	6

Fixation with 4.5 mm cortical locking screws or either 6.5 mm cancellous locking screws



OF 1053.07 L-R	7
OF 1053.07 L-L	7
OF 1053.08 L-R	8
OF 1053.08 L-L	8
OF 1053.09 L-R	9
OF 1053.09 L-R	9
OF 1053.10 L-R	10
OF 1053.10 L-L	10

Lateral Tibial Buttress Plate (R/L)

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

Available in Right and Left version

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 4.5 mm cortical screws or either 6.5 mm cancellous screws



CODE	Holes
OF 1075.05 R	5
OF 1075.05 L	5
OF 1075.06 R	6
OF 1075.06 L	6
OF 1075.07 R	7
OF 1075.07 L	7
OF 1075.08 R	8
OF 1075.08 L	8

OF 1075.09 R	9
OF 1075.09 L	9
OF 1075.10 R	10
OF 1075.10 L	10
OF 1075.11 R	11
OF 1075.11 L	11
OF 1075.12 R	12
OF 1075.12 L	12

Lateral Tibial Buttress Plate - LC (R/L)

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 4.5 mm cortical screws or either 6.5 mm cancellous screws



CODE	Holes
OF 1075.05 R LC	5
OF 1075.05 L LC	5
OF 1075.06 R LC	6
OF 1075.06 L LC	6
OF 1075.07 R LC	7
OF 1075.07 L LC	7
OF 1075.08 R LC	8
OF 1075.08 L LC	8

OF 1075.09 R LC	9
OF 1075.09 L LC	9
OF 1075.10 R LC	10
OF 1075.10 L LC	10
OF 1075.11 R LC	11
OF 1075.11 L LC	11
OF 1075.12 R LC	12
OF 1075.12 L LC	12

Lateral Tibial Buttress Plate (R/L) – Locking

Indication :

A buttress plate applied a forced to the bone which is perpendicular (normal) to the flat surface of the plate

With double bend as buttress for tibial head fractures

Is strengthen (buttress) a weakned area of cortex

Available in Right and Left version

Locking plate

The plate prevents the bone from collapsing during the healing process

Indicated for the osteosynthesis of fractures at tibia at proximal area

Fixation with 4.5 mm cortical locking screws or either 6.5 mm cancellous locking screws



CODE	Holes
OF 1075.05 L-R	5
OF 1075.05 L-L	5
OF 1075.06 L-R	6
OF 1075.06 L-L	6
OF 1075.07 L-R	7
OF 1075.07 L-L	7
OF 1075.08 L-R	8
OF 1075.08 L-L	8

OF 1075.09 L-R	9
OF 1075.09 L-L	9
OF 1075.10 L-R	10
OF 1075.10 L-L	10
OF 1075.11 L-R	11
OF 1075.11 L-L	11
OF 1075.12 L-R	12
OF 1075.12 L-L	12

Multi Condylar Plate (Condylar Buttress) (R/L)

Indication :

As buttress for femoral condyles with multiple fragment fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate

Available in Right and Left version



CODE	Holes
OF 1058.04 R	4
OF 1058.04 L	4
OF 1058.05 R	5

OF 1058.05 L	5
OF 1058.06 R	6
OF 1058.06 L	6
OF 1058.07 R	7



OF 1058.07 L	7
OF 1058.08 R	8
OF 1058.08 L	8
OF 1058.09 R	9
OF 1058.09 L	9
OF 1058.10 R	10
OF 1058.10 L	10
OF 1058.11 R	11
OF 1058.11 L	11

OF 1058.12 R	12
OF 1058.12 L	12
OF 1058.13 R	13
OF 1058.13 L	13
OF 1058.14 R	14
OF 1058.14 L	14
OF 1058.15 R	15
OF 1058.15 L	15

Multi Condylar Plate (Condylar Buttress) (R/L) – LC

Indication :

As buttress for femoral condyles with multiple fragment fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate

Available in Right and Left version



CODE	Holes
OF 1058.04 R	4
OF 1058.04 L	4
OF 1058.05 R	5
OF 1058.05 L	5
OF 1058.06 R	6
OF 1058.06 L	6
OF 1058.07 R	7
OF 1058.07 L	7
OF 1058.08 R	8
OF 1058.08 L	8
OF 1058.09 R	9
OF 1058.09 L	9

OF 1058.10 R	10
OF 1058.10 L	10
OF 1058.11 R	11
OF 1058.11 L	11
OF 1058.12 R	12
OF 1058.12 L	12
OF 1058.13 R	13
OF 1058.13 L	13
OF 1058.14 R	14
OF 1058.14 L	14
OF 1058.15 R	15
OF 1058.15 L	15



Multi Condylar Plate (Condylar Buttress) (R/L) – Locking

Indication :

As buttress for femoral condyles with multiple fragment fractures

Is strengthen (buttress) a weakned area of cortex

The plate prevents the bone from collapsing during the healing process

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

Fixation with 4.5 mm cortical locking screws, 6.5 mm cancellous locking screws are used in the head part of plate

Available in Right and Left version



CODE	Holes
OF 1058.04 L-R	4
OF 1058.04 L-L	4
OF 1058.05 L-R	5
OF 1058.05 L-L	5

OF 1058.06 L-R	6
OF 1058.06 L-L	6
OF 1058.07 L-R	7
OF 1058.07 L-L	7
OF 1058.08 L-R	8
OF 1058.08 L-L	8
OF 1058.09 L-R	9
OF 1058.09 L-L	9
OF 1058.10 L-R	10
OF 1058.10 L-L	10
OF 1058.11 L-R	11
OF 1058.11 L-L	11
OF 1058.12 L-R	12
OF 1058.12 L-L	12
OF 1058.13 -LR	13
OF 1058.13 L-L	13
OF 1058.14 L-R	14
OF 1058.14 L-L	14
OF 1058.15 L-R	15
OF 1058.15 L-L	15

Trochanter Ear Plate

Indication :

Used for proximal femur fractures

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1072.04 A	4
OF 1072.05 A	5
OF 1072.06 A	6

OF 1072.07 A	7
OF 1072.08 A	8
OF 1072.09 A	9
OF 1072.10 A	10



OF 1072.11 A	11
OF 1072.12 A	12
OF 1072.13 A	13

OF 1072.14 A	14
OF 1072.15 A	15

Trochanter Ear Plate – LC

Indication :

Used for proximal femur fractures

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1072.04 A LC	4
OF 1072.05 A LC	5
OF 1072.06 A LC	6
OF 1072.07 A LC	7
OF 1072.08 A LC	8
OF 1072.09 A LC	9

OF 1072.10 A LC	10
OF 1072.11 A LC	11
OF 1072.12 A LC	12
OF 1072.13 A LC	13
OF 1072.14 A LC	14
OF 1072.15 A LC	15

Trochanter Ear Plate – Locking

Indication :

Used for proximal femur fractures

Fixation with 4.5 mm cortical locking screws, 6.5 mm cancellous locking screws are used in the head part of plate



CODE	Holes
OF 1072.04 AL	4
OF 1072.05 AL	5
OF 1072.06 AL	6
OF 1072.07 AL	7
OF 1072.08 AL	8

OF 1072.09 AL	9
OF 1072.10 AL	10
OF 1072.11 AL	11
OF 1072.12 AL	12
OF 1072.13 AL	13
OF 1072.14 AL	14
OF 1072.15 AL	15



Trochanter Straight Plate

Indication :

Used for inter sub – trochanteric fractures

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1072.04 B	4
OF 1072.05 B	5
OF 1072.06 B	6
OF 1072.07 B	7
OF 1072.08 B	8
OF 1072.09 B	9
OF 1072.10 B	10
OF 1072.11 B	11
OF 1072.12 B	12
OF 1072.13 B	13
OF 1072.14 B	14
OF 1072.15 B	15

Trochanter Straight Plate – LC

Indication :

Used for inter sub – trochanteric fractures

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1072.04 B LC	4
OF 1072.05 B LC	5
OF 1072.06 B LC	6
OF 1072.07 B LC	7
OF 1072.08 B LC	8
OF 1072.09 B LC	9
OF 1072.10 B LC	10
OF 1072.11 B LC	11
OF 1072.12 B LC	12
OF 1072.13 B LC	13
OF 1072.14 B LC	14
OF 1072.15 B LC	15



Trochanter Straight Plate – Locking

Indication :

Used for inter sub – trochanteric fractures

Fixation with 4.5 mm cortical locking screws, 6.5 mm cancellous locking screws are used in the head part of plate



CODE	Holes
OF 1072.04 BL	4
OF 1072.05 BL	5
OF 1072.06 BL	6
OF 1072.07 BL	7
OF 1072.08 BL	8
OF 1072.09 BL	9
OF 1072.10 BL	10
OF 1072.11 BL	11
OF 1072.12 BL	12
OF 1072.13 BL	13
OF 1072.14 BL	14
OF 1072.15 BL	15

Cobra Plate

Indication :

For arthrodesis of the hip shaft with dynamic compression holes

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1064.04	4
OF 1064.05	5
OF 1064.06	6
OF 1064.07	7
OF 1064.08	8
OF 1064.09	9
OF 1064.10	10
OF 1064.11	11
OF 1064.12	12
OF 1064.13	13
OF 1064.14	14
OF 1064.15	15



Cobra Plate – LC

Indication :

For arthrodesis of the hip shaft with dynamic compression holes

Greatly reduced the area of the plate bone contact

Less compromised the capillary network of periosteum

Reduce the porotic changes underneath the plate

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1064.04 LC	4
OF 1064.05 LC	5
OF 1064.06 LC	6
OF 1064.07 LC	7
OF 1064.08 LC	8
OF 1064.09 LC	9
OF 1064.10 LC	10
OF 1064.11 LC	11
OF 1064.12 LC	12
OF 1064.13 LC	13
OF 1064.14 LC	14
OF 1064.15 LC	15

Cobra Plate – Locking

Indication :

For arthrodesis of the hip shaft with dynamic compression holes

Fixation with 4.5 mm cortical locking screws, 6.5 mm cancellous locking screws are used in the head part of plate



CODE	Holes
OF 1064.04 L	4
OF 1064.05 L	5
OF 1064.06 L	6
OF 1064.07 L	7
OF 1064.08 L	8
OF 1064.09 L	9

OF 1064.10 L	10
OF 1064.11 L	11
OF 1064.12 L	12
OF 1064.13 L	13
OF 1064.14 L	14
OF 1064.15 L	15



Distal Femoral Plate (R/L)

Indication :

Indicated for for
supracondylar, intra – articular
and extra – articular condylar
fractures

Fractures in osteopenic bone,
non-unions, mal-unions

Fixation with 4.5 mm cortical
screws, 6.5 mm cancellous
screws are used in the head
part of plate



CODE	Holes
OF 1074.04 R	4
OF 1074.04 L	4
OF 1074.05 R	5
OF 1074.05 L	5
OF 1074.06 R	6
OF 1074.06 L	6
OF 1074.07 R	7
OF 1074.07 L	7
OF 1074.08 R	8
OF 1074.08 L	8
OF 1074.09 R	9
OF 1074.09 L	9

OF 1074.10 R	10
OF 1074.10 L	10
OF 1074.11 R	11
OF 1074.11 L	11
OF 1074.12 R	12
OF 1074.12 L	12
OF 1074.13 R	13
OF 1074.13 L	13
OF 1074.14 R	14
OF 1074.14 L	14
OF 1074.15 R	15
OF 1074.15 L	15

Distal Femoral Plate (R/L) – L

Indication :

Indicated for for
supracondylar, intra – articular
and extra – articular condylar
fractures

Greatly reduced the area of
the plate bone contact

Less compromised the
capillary network of
periosteum

Reduce the porotic changes
underneath the plate

Fractures in osteopenic bone,
non-unions, mal-unions

Fixation with 4.5 mm cortical
screws, 6.5 mm cancellous
screws are used in the head
part of plate



CODE	Holes
OF 1074.04 R LC	4
OF 1074.04 L LC	4
OF 1074.05 R LC	5
OF 1074.05 L LC	5
OF 1074.06 R LC	6
OF 1074.06 L LC	6
OF 1074.07 R LC	7
OF 1074.07 L LC	7

OF 1074.08 R LC	8
OF 1074.08 L LC	8
OF 1074.09 R LC	9
OF 1074.09 L LC	9
OF 1074.10 R LC	10
OF 1074.10 L LC	10
OF 1074.11 R LC	11
OF 1074.11 L LC	11
OF 1074.12 R LC	12

OF 1074.12 L LC	12
OF 1074.13 R LC	13
OF 1074.13 L LC	13
OF 1074.14 R LC	14
OF 1074.14 L LC	14
OF 1074.15 R LC	15
OF 1074.15 L LC	15

Distal Femoral Plate (R/L) – Locking

Indication :

Indicated for for
supracondylar, intra – articular
and extra – articular condylar
fractures

Fractures in osteopenic bone,
non-unions, mal-unions

Fixation with 4.5 mm cortical
screws, 6.5 mm cancellous
screws are used in the head
part of plate





CODE	Holes
OF 1074.04 LR	4
OF 1074.04 LL	4
OF 1074.05 LR	5
OF 1074.05 LL	5
OF 1074.06 LR	6
OF 1074.06 LL	6
OF 1074.07 LR	7
OF 1074.07 LL	7

OF 1074.08 LR	8
OF 1074.08 LL	8
OF 1074.09 LR	9
OF 1074.09 LL	9
OF 1074.10 LR	10
OF 1074.10 LL	10
OF 1074.11 LR	11
OF 1074.11 LL	11
OF 1074.12 LR	12

OF 1074.12 LL	12
OF 1074.13 LR	13
OF 1074.13 LL	13
OF 1074.14 LR	14
OF 1074.14 LL	14
OF 1074.15 LR	15
OF 1074.15 LL	15

Proximal Lateral Femoral Plate (R/L)

Indication :

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1082.04 R	4
OF 1082.04 L	4
OF 1082.05 R	5
OF 1082.05 L	5
OF 1082.06 R	6
OF 1082.06 L	6
OF 1082.07 R	7
OF 1082.07 L	7

OF 1082.08 R	8
OF 1082.08 L	8
OF 1082.09 R	9
OF 1082.09 L	9
OF 1082.10 R	10
OF 1082.10 L	10
OF 1082.11 R	11
OF 1082.11 L	11
OF 1082.12 R	12

OF 1082.12 L	12
OF 1082.13 R	13
OF 1082.13 L	13
OF 1082.14 R	14
OF 1082.14 L	14
OF 1082.15 R	15
OF 1082.15 L	15



Proximal Lateral Femoral Plate (R/L) – LC

Indication :

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1082.04 R LC	4
OF 1082.04 L LC	4
OF 1082.05 R LC	5
OF 1082.05 L LC	5
OF 1082.06 R LC	6
OF 1082.06 L LC	6
OF 1082.07 R LC	7
OF 1082.07 L LC	7

OF 1082.08 R LC	8
OF 1082.08 L LC	8
OF 1082.09 R LC	9
OF 1082.09 L LC	9
OF 1082.10 R LC	10
OF 1082.10 L LC	10
OF 1082.11 R LC	11
OF 1082.11 L LC	11
OF 1082.12 R LC	12

OF 1082.12 L LC	12
OF 1082.13 R LC	13
OF 1082.13 L LC	13
OF 1082.14 R LC	14
OF 1082.14 L LC	14
OF 1082.15 R LC	15
OF 1082.15 L LC	15

Proximal Lateral Femoral Plate (R/L) – Locking

Indication :

Fixation with 4.5 mm cortical locking screws, 6.5 mm cancellous locking screws are used in the head part of plate



CODE	Holes
OF 1082.04 LR	4
OF 1082.04 LL	4
OF 1082.05 LR	5
OF 1082.05 LL	5
OF 1082.06 LR	6
OF 1082.06 LL	6
OF 1082.07 LR	7
OF 1082.07 LL	7

OF 1082.08 LR	8
OF 1082.08 LL	8
OF 1082.09 LR	9
OF 1082.09 LL	9
OF 1082.10 LR	10
OF 1082.10 LL	10
OF 1082.11 LR	11
OF 1082.11 LL	11
OF 1082.12 LR	12

OF 1082.12 LL	12
OF 1082.13 LR	13
OF 1082.13 LL	13
OF 1082.14 LR	14
OF 1082.14 LL	14
OF 1082.15 LR	15
OF 1082.15 LL	15

Proximal Lateral Tibia – Golf Plate (R/L)

Indication :

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1084.04 R	4
OF 1084.04 L	4
OF 1084.05 R	5
OF 1084.05 L	5
OF 1084.06 R	6
OF 1084.06 L	6
OF 1084.07 R	7
OF 1084.07 L	7

OF 1084.08 R	8
OF 1084.08 L	8
OF 1084.09 R	9
OF 1084.09 L	9
OF 1084.10 R	10
OF 1084.10 L	10
OF 1084.11 R	11
OF 1084.11 L	11
OF 1084.12 R	12

OF 1084.12 L	12
OF 1084.13 R	13
OF 1084.13 L	13
OF 1084.14 R	14
OF 1084.14 L	14
OF 1084.15 R	15
OF 1084.15 L	15

Proximal Lateral Tibia – Golf Plate (R/L) – LC

Indication :

Fixation with 4.5 mm cortical screws, 6.5 mm cancellous screws are used in the head part of plate



CODE	Holes
OF 1084.04 R LC	4
OF 1084.04 L LC	4
OF 1084.05 R LC	5
OF 1084.05 L LC	5
OF 1084.06 R LC	6
OF 1084.06 L LC	6
OF 1084.07 R LC	7
OF 1084.07 L LC	7

OF 1084.08 R LC	8
OF 1084.08 L LC	8
OF 1084.09 R LC	9
OF 1084.09 L LC	9
OF 1084.10 R LC	10
OF 1084.10 L LC	10
OF 1084.11 R LC	11
OF 1084.11 L LC	11
OF 1084.12 R LC	12

OF 1084.12 L LC	12
OF 1084.13 R LC	13
OF 1084.13 L LC	13
OF 1084.14 R LC	14
OF 1084.14 L LC	14
OF 1084.15 R LC	15
OF 1084.15 L LC	15

Proximal Lateral Tibia – Golf Plate (R/L) – Locking

Indication :

Fixation with 4.5 mm cortical locking screws, 6.5 mm cancellous locking screws are used in the head part of plate



CODE	Holes
OF 1084.04 LR	4
OF 1084.04 LL	4
OF 1084.05 LR	5
OF 1084.05 LL	5
OF 1084.06 LR	6
OF 1084.06 LL	6
OF 1084.07 LR	7
OF 1084.07 LL	7

OF 1084.08 LR	8
OF 1084.08 LL	8
OF 1084.09 LR	9
OF 1084.09 LL	9
OF 1084.10 LR	10
OF 1084.10 LL	10
OF 1084.11 LR	11
OF 1084.11 LL	11
OF 1084.12 LR	12

OF 1084.12 LL	12
OF 1084.13 LR	13
OF 1084.13 LL	13
OF 1084.14 LR	14
OF 1084.14 LL	14
OF 1084.15 LR	15
OF 1084.15 LL	15

4.5 mm Reconstruction Plate – Straight

Indication :

Accurate contouring in any plane

Used for calcaneus and pelvic acetabulum fractures

Deep notches between holes



Fixation with 4.5 mm cortical screw

CODE	Holes
OF 1054.04	4
OF 1054.05	5
OF 1054.06	6
OF 1054.07	7
OF 1054.08	8
OF 1054.09	9

OF 1054.10	10
OF 1054.11	11
OF 1054.12	12
OF 1054.13	13
OF 1054.14	14
OF 1054.15	15



4.5 mm Reconstruction Plate – Straight – Locking

Indication :

Accurate contouring in any plane

Used for calcaneus and pelvic acetabulum fractures

Deep notches between holes



Fixation with 4.5 mm cortical locking screw

CODE	Holes
OF 1054.04 L	4
OF 1054.05 L	5
OF 1054.06 L	6
OF 1054.07 L	7
OF 1054.08 L	8
OF 1054.09 L	9
OF 1054.10 L	10
OF 1054.11 L	11
OF 1054.12 L	12
OF 1054.13 L	13
OF 1054.14 L	14
OF 1054.15 L	15

4.5 mm Reconstruction Plate – Curved

Indication :

Used for pelvic and acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 4.5 mm cortical screw



Plate radius : 100

CODE	Holes
OF 1065.04	4
OF 1065.05	5
OF 1065.06	6
OF 1065.07	7
OF 1065.08	8
OF 1065.09	9

OF 1065.10	10
OF 1065.11	11
OF 1065.12	12
OF 1065.13	13
OF 1065.14	14
OF 1065.15	15



4.5 mm Reconstruction Plate – Curved – Locking

Indication :

Used for pelvic and acetabulum fractures

Deep notches between holes

Accurate contouring in any plane

Fixation with 4.5 mm cortical locking screw

Plate radius : 100



CODE	Holes
OF 1065.04 L	4
OF 1065.05 L	5
OF 1065.06 L	6
OF 1065.07 L	7
OF 1065.08 L	8
OF 1065.09 L	9
OF 1065.10 L	10
OF 1065.11 L	11
OF 1065.12 L	12
OF 1065.13 L	13
OF 1065.14 L	14
OF 1065.15 L	15

4.5 mm Sherman Plate – Straight

Indication :

Used for pelvic acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 4.5 mm cortical screw

Plate radius : 108



CODE	Holes
OF 1055.04	4
OF 1055.05	5
OF 1055.06	6
OF 1055.07	7
OF 1055.08	8
OF 1055.09	9

OF 1055.10	10
OF 1055.11	11
OF 1055.12	12
OF 1055.13	13
OF 1055.14	14
OF 1055.15	15



4.5 mm Sherman Plate – Straight – Locking

Indication :

Used for pelvic acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 4.5 mm cortical locking screw

Plate radius : 108



CODE	Holes
OF 1055.04 L	4
OF 1055.05 L	5
OF 1055.06 L	6
OF 1055.07 L	7
OF 1055.08 L	8
OF 1055.09 L	9
OF 1055.10 L	10
OF 1055.11 L	11
OF 1055.12 L	12
OF 1055.13 L	13
OF 1055.14 L	14
OF 1055.15 L	15

4.5 mm Sherman Plate – Curved

Indication :

Used for pelvic acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 4.5 mm cortical screw

Plate radius : 108



CODE	Holes
OF 1056.04	4
OF 1056.05	5
OF 1056.06	6
OF 1056.07	7
OF 1056.08	8
OF 1056.09	9

OF 1056.10	10
OF 1056.11	11
OF 1056.12	12
OF 1056.13	13
OF 1056.14	14
OF 1056.15	15



4.5 mm Sherman Plate – Curved – Locking

Indication :

Used for pelvic acetabular surgery

Deep notches between holes

Accurate contouring in any plane

Fixation with 4.5 mm cortical screw



CODE	Holes
OF 1056.04 L	4
OF 1056.05 L	5
OF 1056.06 L	6
OF 1056.07 L	7
OF 1056.08 L	8
OF 1056.09 L	9
OF 1056.10 L	10
OF 1056.11 L	11
OF 1056.12 L	12
OF 1056.13 L	13
OF 1056.14 L	14
OF 1056.15 L	15

Spoon Plate (V Plate)

Indication :

Used for anterior crest of the proximal tibia

Shaft with V profile for the fixation on the frontal edge of the tibia

Fixation with 4.5 mm cortical screw and 6.5 mm cancellous screw in the head part of plate



CODE	Holes
OF 1057.05	5
OF 1057.06	6

OF 1057.07	7
OF 1057.08	8
OF 1057.09	9

OF 1057.10	10
OF 1057.11	11
OF 1057.12	12



4.5 mm Hooked Plate DCP

Indication :

Used for proximal humerus, olecranon

Fixation with 4.5 mm cortical screw



CODE	Holes
OF 10671.03	3
OF 10671.04	4
OF 10671.05	5
OF 10671.06	6
OF 10671.07	7
OF 10671.08	8
OF 10671.09	9
OF 10671.10	10

DHS Plate (Dynamic Hip System)

Indication :

Used for petrochanteric and intertrochanteric fracture (proximal femur)

Fixation with 4.5 mm cortical screw and lag screw

Plate angle : 130°

Barrel length : 38 mm

Outer barrel dia : 12.6 mm



CODE	Holes
OF 4020.03	3
OF 4020.04	4
OF 4020.05	5
OF 4020.06	6
OF 4020.07	7

OF 4020.08	8
OF 4020.09	9
OF 4020.10	10
OF 4020.11	11
OF 4020.12	12



DHS Plate (Dynamic Hip System) – LC

Indication :

Used for petrochanteric and intertrochanteric fracture (proximal femur)

Fixation with 4.5 mm cortical screw and lag screw

Plate angle : 130°

Barrel length : 38 mm

Outer barrel dia : 12.6 mm



CODE	Holes
OF 4020.03 LC	3
OF 4020.04 LC	4
OF 4020.05 LC	5
OF 4020.06 LC	6
OF 4020.07 LC	7
OF 4020.08 LC	8
OF 4020.09 LC	9
OF 4020.10 LC	10
OF 4020.11 LC	11
OF 4020.12 LC	12

DCS Plate (Dynamic Condyle System)

Indication :

For gentle tretment of condylar fracture

Fixation with 4.5 mm cortical screw and lag screw

Plate angle : 95°

Barrel length : 34 mm

Outer barrel dia : 12.6 mm



CODE	Holes
OF 4040.03	3
OF 4040.04	4
OF 4040.05	5

OF 4040.06	6
OF 4040.07	7
OF 4040.08	8
OF 4040.09	9

OF 4040.10	10
OF 4040.11	11
OF 4040.12	12



DCS Plate (Dynamic Condyle System) – LC

Indication :

Used for petrochanteric and intertrochanteric fracture (proximal femur)

Fixation with 4.5 mm cortical screw and lag screw

Plate angle : 130°

Barrel length : 38 mm

Outer barrel dia : 12.6 mm



CODE	Holes
OF 4040.03 LC	3
OF 4040.04 LC	4
OF 4040.05 LC	5
OF 4040.06 LC	6
OF 4040.07 LC	7
OF 4040.08 LC	8
OF 4040.09 LC	9
OF 4040.10 LC	10
OF 4040.11 LC	11
OF 4040.12 LC	12



4. ANGULATED BLADE PLATE

High Tibia Osteotomy Plate (R/L)

Indication :

Used for tibia head, distal femur for small fragments

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Overall height : 94 mm

Plate angle : 90°

Available in Right and Left version

Blade profile H



CODE						Holes
Blade 50	Blade 55	Blade 60	Blade 65	Blade 70	Blade 75	
OF 4050.04	OF 4055.04	OF 4060.04	OF 4065.04	OF 4070.04	OF 4075.04	4
OF 4050.05	OF 4055.05	OF 4060.05	OF 4065.05	OF 4070.05	OF 4075.05	5
OF 4050.06	OF 4055.06	OF 4060.06	OF 4065.06	OF 4070.06	OF 4075.06	6
OF 4050.07	OF 4055.07	OF 4060.07	OF 4065.07	OF 4070.07	OF 4075.07	7
OF 4050.08	OF 4055.08	OF 4060.08	OF 4065.08	OF 4070.08	OF 4075.08	8
OF 4050.09	OF 4055.09	OF 4060.09	OF 4065.09	OF 4070.09	OF 4075.09	9
OF 4050.10	OF 4055.10	OF 4060.10	OF 4065.10	OF 4070.10	OF 4075.10	10

High Tibia Osteotomy Plate (R/L) – LC

Indication :

Used for tibia head, distal femur for small fragments

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Plate angle : 90°

Available in Right and Left version

Blade profile H with low contact plate



:

CODE						Holes
Blade 50	Blade 55	Blade 60	Blade 65	Blade 70	Blade 75	
OF 4050.04 LC	OF 4055.04 LC	OF 4060.04 LC	OF 4065.04 LC	OF 4070.04 LC	OF 4075.04 LC	4
OF 4050.05 LC	OF 4055.05 LC	OF 4060.05 LC	OF 4065.05 LC	OF 4070.05 LC	OF 4075.05 LC	5
OF 4050.06 LC	OF 4055.06 LC	OF 4060.06 LC	OF 4065.06 LC	OF 4070.06 LC	OF 4075.06 LC	6
OF 4050.07 LC	OF 4055.07 LC	OF 4060.07 LC	OF 4065.07 LC	OF 4070.07 LC	OF 4075.07 LC	7
OF 4050.08 LC	OF 4055.08 LC	OF 4060.08 LC	OF 4065.08 LC	OF 4070.08 LC	OF 4075.08 LC	8
OF 4050.09 LC	OF 4055.09 LC	OF 4060.09 LC	OF 4065.09 LC	OF 4070.09 LC	OF 4075.09 LC	9
OF 4050.10 LC	OF 4055.10 LC	OF 4060.10 LC	OF 4065.10 LC	OF 4070.10 LC	OF 4075.10 LC	10

Petrochanteric Plate

Indication :

Used for tibia head, distal femur for small fragments

Fixation with 4.5 mm cortical screw or either 6.5 mm cancellous screw

Plate angle : 135°

Blade profile H



CODE						Holes
Blade 50	Blade 55	Blade 60	Blade 65	Blade 70	Blade 75	
OF 5050.04	OF 5055.04	OF 5060.04	OF 5065.04	OF 5070.04	OF 5075.04	4
OF 5050.05	OF 5055.05	OF 5060.05	OF 5065.05	OF 5070.05	OF 5075.05	5
OF 5050.06	OF 5055.06	OF 5060.06	OF 5065.06	OF 5070.06	OF 5075.06	6
OF 5050.07	OF 5055.07	OF 5060.07	OF 5065.07	OF 5070.07	OF 5075.07	7
OF 5050.08	OF 5055.08	OF 5060.08	OF 5065.08	OF 5070.08	OF 5075.08	8



OF 5050.09	OF 5055.09	OF 5060.09	OF 5065.09	OF 5070.09	OF 5075.09	9
OF 5050.10	OF 5055.10	OF 5060.10	OF 5065.10	OF 5070.10	OF 5075.10	10
OF 5050.11	OF 5055.11	OF 5060.11	OF 5065.11	OF 5070.11	OF 5075.11	11
OF 5050.12	OF 5055.12	OF 5060.12	OF 5065.12	OF 5070.12	OF 5075.12	12
OF 5050.13	OF 5055.13	OF 5060.13	OF 5065.13	OF 5070.13	OF 5075.13	13
OF 5050.14	OF 5055.14	OF 5060.14	OF 5065.14	OF 5070.14	OF 5075.14	14
OF 5050.15	OF 5055.15	OF 5060.15	OF 5065.15	OF 5070.15	OF 5075.15	15
OF 5050.16	OF 5055.16	OF 5060.16	OF 5065.16	OF 5070.16	OF 5075.16	16
OF 5050.17	OF 5055.17	OF 5060.17	OF 5065.17	OF 5070.17	OF 5075.17	17
OF 5050.18	OF 5055.18	OF 5060.18	OF 5065.18	OF 5070.18	OF 5075.18	18
OF 5050.19	OF 5055.19	OF 5060.19	OF 5065.19	OF 5070.19	OF 5075.19	19
OF 5050.20	OF 5055.20	OF 5060.20	OF 5065.20	OF 5070.20	OF 5075.20	20

Petrochanteric Plate – LC

Indication :

Used for tibia head, distal
femur for small fragments

Fixation with 4.5 mm cortical
screw or either 6.5 mm
cancellous screw

Plate angle : 135°

Blade profile H



CODE						Holes
Blade 50	Blade 55	Blade 60	Blade 65	Blade 70	Blade 75	
OF 5050.04	OF 5055.04	OF 5060.04	OF 5065.04	OF 5070.04	OF 5075.04	4
OF 5050.05	OF 5055.05	OF 5060.05	OF 5065.05	OF 5070.05	OF 5075.05	5
OF 5050.06	OF 5055.06	OF 5060.06	OF 5065.06	OF 5070.06	OF 5075.06	6
OF 5050.07	OF 5055.07	OF 5060.07	OF 5065.07	OF 5070.07	OF 5075.07	7
OF 5050.08	OF 5055.08	OF 5060.08	OF 5065.08	OF 5070.08	OF 5075.08	8
OF 5050.09	OF 5055.09	OF 5060.09	OF 5065.09	OF 5070.09	OF 5075.09	9
OF 5050.10	OF 5055.10	OF 5060.10	OF 5065.10	OF 5070.10	OF 5075.10	10
OF 5050.11	OF 5055.11	OF 5060.11	OF 5065.11	OF 5070.11	OF 5075.11	11
OF 5050.12	OF 5055.12	OF 5060.12	OF 5065.12	OF 5070.12	OF 5075.12	12
OF 5050.13	OF 5055.13	OF 5060.13	OF 5065.13	OF 5070.13	OF 5075.13	13
OF 5050.14	OF 5055.14	OF 5060.14	OF 5065.14	OF 5070.14	OF 5075.14	14
OF 5050.15	OF 5055.15	OF 5060.15	OF 5065.15	OF 5070.15	OF 5075.15	15
OF 5050.16	OF 5055.16	OF 5060.16	OF 5065.16	OF 5070.16	OF 5075.16	16
OF 5050.17	OF 5055.17	OF 5060.17	OF 5065.17	OF 5070.17	OF 5075.17	17
OF 5050.18	OF 5055.18	OF 5060.18	OF 5065.18	OF 5070.18	OF 5075.18	18
OF 5050.19	OF 5055.19	OF 5060.19	OF 5065.19	OF 5070.19	OF 5075.19	19
OF 5050.20	OF 5055.20	OF 5060.20	OF 5065.20	OF 5070.20	OF 5075.20	20



Condylar ABP Plate 95°

Indication :

Used fracture of distal and proximal femur, proximal tibia, intertrochanteric valgus osteotomy

Fixation with 4.5 mm cortical screw, the two holes next to the blade are larger to accept 6.5 mm cancellous screw

Plate angle : 95°

Blade profile H



CODE						Holes
Blade 50	Blade 55	Blade 60	Blade 65	Blade 70	Blade 75	
OF 4150.04	OF 4155.04	OF 4160.04	OF 4165.04	OF 4170.04	OF 4175.04	4
OF 4150.05	OF 4155.05	OF 4160.05	OF 4165.05	OF 4170.05	OF 4175.05	5
OF 4150.06	OF 4155.06	OF 4160.06	OF 4165.06	OF 4170.06	OF 4175.06	6
OF 4150.07	OF 4155.07	OF 4160.07	OF 4165.07	OF 4170.07	OF 4175.07	7
OF 4150.08	OF 4155.08	OF 4160.08	OF 4165.08	OF 4170.08	OF 4175.08	8
OF 4150.09	OF 4155.09	OF 4160.09	OF 4165.09	OF 4170.09	OF 4175.09	9
OF 4150.10	OF 4155.10	OF 4160.10	OF 4165.10	OF 4170.10	OF 4175.10	10
OF 4150.11	OF 4155.11	OF 4160.11	OF 4165.11	OF 4170.11	OF 4175.11	11
OF 4150.12	OF 4155.12	OF 4160.12	OF 4165.12	OF 4170.12	OF 4175.12	12
OF 4150.13	OF 4155.13	OF 4160.13	OF 4165.13	OF 4170.13	OF 4175.13	13
OF 4150.14	OF 4155.14	OF 4160.14	OF 4165.14	OF 4170.14	OF 4175.14	14
OF 4150.15	OF 4155.15	OF 4160.15	OF 4165.15	OF 4170.15	OF 4175.15	15
OF 4150.16	OF 4155.16	OF 4160.16	OF 4165.16	OF 4170.16	OF 4175.16	16
OF 4150.17	OF 4155.17	OF 4160.17	OF 4165.17	OF 4170.17	OF 4175.17	17
OF 4150.18	OF 4155.18	OF 4160.18	OF 4165.18	OF 4170.18	OF 4175.18	18
OF 4150.19	OF 4155.19	OF 4160.19	OF 4165.19	OF 4170.19	OF 4175.19	19
OF 4150.20	OF 4155.20	OF 4160.20	OF 4165.20	OF 4170.20	OF 4175.20	20

Condylar ABP Plate 95° – LC

Indication :

Used fracture of distal and proximal femur, proximal tibia, intertrochanteric valgus osteotomy

Fixation with 4.5 mm cortical screw, the two holes next to the blade are larger to accept 6.5 mm cancellous screw

Plate angle : 95°

Blade profile H





CODE						Holes
Blade 50	Blade 55	Blade 60	Blade 65	Blade 70	Blade 75	
OF 4150.04	OF 4155.04	OF 4160.04	OF 4165.04	OF 4170.04	OF 4175.04	4
OF 4150.05	OF 4155.05	OF 4160.05	OF 4165.05	OF 4170.05	OF 4175.05	5
OF 4150.06	OF 4155.06	OF 4160.06	OF 4165.06	OF 4170.06	OF 4175.06	6
OF 4150.07	OF 4155.07	OF 4160.07	OF 4165.07	OF 4170.07	OF 4175.07	7
OF 4150.08	OF 4155.08	OF 4160.08	OF 4165.08	OF 4170.08	OF 4175.08	8
OF 4150.09	OF 4155.09	OF 4160.09	OF 4165.09	OF 4170.09	OF 4175.09	9
OF 4150.10	OF 4155.10	OF 4160.10	OF 4165.10	OF 4170.10	OF 4175.10	10
OF 4150.11	OF 4155.11	OF 4160.11	OF 4165.11	OF 4170.11	OF 4175.11	11
OF 4150.12	OF 4155.12	OF 4160.12	OF 4165.12	OF 4170.12	OF 4175.12	12
OF 4150.13	OF 4155.13	OF 4160.13	OF 4165.13	OF 4170.13	OF 4175.13	13
OF 4150.14	OF 4155.14	OF 4160.14	OF 4165.14	OF 4170.14	OF 4175.14	14
OF 4150.15	OF 4155.15	OF 4160.15	OF 4165.15	OF 4170.15	OF 4175.15	15
OF 4150.16	OF 4155.16	OF 4160.16	OF 4165.16	OF 4170.16	OF 4175.16	16
OF 4150.17	OF 4155.17	OF 4160.17	OF 4165.17	OF 4170.17	OF 4175.17	17
OF 4150.18	OF 4155.18	OF 4160.18	OF 4165.18	OF 4170.18	OF 4175.18	18
OF 4150.19	OF 4155.19	OF 4160.19	OF 4165.19	OF 4170.19	OF 4175.19	19
OF 4150.20	OF 4155.20	OF 4160.20	OF 4165.20	OF 4170.20	OF 4175.20	20

Small ABP Plate 80° /90°/ 100°

Indication :

Used fracture of distal and proximal femur, proximal tibia, intertrochanteric valgus osteotomy for children

Fixation with 4.5 mm cortical screw

Plate angle 80° /90° /100°

Blade profile H



CODE						Holes
Blade 35	Blade 40	Blade 45	Blade 50	Blade 55	Blade 60	
OF 4935.04	OF 4940.04	OF 4945.04	OF 4950.04	OF 4955.04	OF 4960.04	4
OF 4935.05	OF 4940.05	OF 4945.05	OF 4950.05	OF 4955.05	OF 4960.05	5
OF 4935.06	OF 4940.06	OF 4945.06	OF 4950.06	OF 4955.06	OF 4960.06	6
OF 4935.07	OF 4940.07	OF 4945.07	OF 4950.07	OF 4955.07	OF 4960.07	7
OF 4935.08	OF 4940.08	OF 4945.08	OF 4950.08	OF 4955.08	OF 4960.08	8
OF 4935.09	OF 4940.09	OF 4945.09	OF 4950.09	OF 4955.09	OF 4960.09	9
OF 4935.10	OF 4940.10	OF 4945.10	OF 4950.10	OF 4955.10	OF 4960.10	10
OF 4935.11	OF 4940.11	OF 4945.11	OF 4950.11	OF 4955.11	OF 4960.11	11
OF 4935.12	OF 4940.12	OF 4945.12	OF 4950.12	OF 4955.12	OF 4960.12	12



Small ABP Plate 80° /90°/ 100° - LC

Indication :

Used fracture of distal and proximal femur, proximal tibia, intertrochanteric valgus osteotomy for children

Fixation with 4.5 mm cortical screw

Plate angle 80° /90° /100°

Blade profile H with low contact plate



CODE						Holes
Blade 35	Blade 40	Blade 45	Blade 50	Blade 55	Blade 60	
OF 4935.04 LC	OF 4940.04 LC	OF 4945.04 LC	OF 4950.04 LC	OF 4955.04 LC	OF 4960.04 LC	4
OF 4935.05 LC	OF 4940.05 LC	OF 4945.05 LC	OF 4950.05 LC	OF 4955.05 LC	OF 4960.05 LC	5
OF 4935.06 LC	OF 4940.06 LC	OF 4945.06 LC	OF 4950.06 LC	OF 4955.06 LC	OF 4960.06 LC	6
OF 4935.07 LC	OF 4940.07 LC	OF 4945.07 LC	OF 4950.07 LC	OF 4955.07 LC	OF 4960.07 LC	7
OF 4935.08 LC	OF 4940.08 LC	OF 4945.08 LC	OF 4950.08 LC	OF 4955.08 LC	OF 4960.08 LC	8
OF 4935.09 LC	OF 4940.09 LC	OF 4945.09 LC	OF 4950.09 LC	OF 4955.09 LC	OF 4960.09 LC	9
OF 4935.10 LC	OF 4940.10 LC	OF 4945.10 LC	OF 4950.10 LC	OF 4955.10 LC	OF 4960.10 LC	10
OF 4935.11 LC	OF 4940.11 LC	OF 4945.11 LC	OF 4950.11 LC	OF 4955.11 LC	OF 4960.11 LC	11
OF 4935.12 LC	OF 4940.12 LC	OF 4945.12 LC	OF 4950.12 LC	OF 4955.12 LC	OF 4960.12 LC	12



5. SPINE PLATE

Cervical Vertebra Plate

Indication :

Used for fractures and spondylodeses spine surgery (vertebral column)

Fixation with 3.5 mm cortical screw



CODE	Holes
OF 1060.05	5
OF 1060.08	8
OF 1060.11	11
OF 1060.14	14
OF 1060.16	16

Stefee Plate

Indication :

Used for fractures and spondylodeses spine surgery (vertebral column)

Fixation with stefee screw



CODE	Holes
OF 1160-1	1
OF 1160-1.5	1.5
OF 1160-2	2
OF 1160-2.5	2.5
OF 1160-3	3
OF 1160-3.5	3.5
OF 1160-4	4




6. SCREW

3.5 mm CORTICAL SCREW

Type of bone screw : small screws for small bone and fragment

Fully threaded

Spherical head with small hexagonal socket



CODE	Length (mm)
OF 510 08	8
OF 510 10	10
OF 510 12	12
OF 510 14	14
OF 510 16	16
OF 510 18	18
OF 510 20	20
OF 510 22	22
OF 510 24	24
OF 510 26	26
OF 510 28	28
OF 510 30	30
OF 510 32	32


OF 510 34	34
OF 510 36	36
OF 510 38	38
OF 510 40	40
OF 510 42	42
OF 510 44	44
OF 510 46	46
OF 510 48	48
OF 510 50	50
OF 510 52	52
OF 510 54	54
OF 510 56	56
OF 510 58	58
OF 510 60	60

4.5 mm CORTICAL SCREW

Type of bone screw : small screws for large bone and fragment

Fully threaded

Spherical head with small hexagonal socket





CODE	Length (mm)
OF 500 10	10
OF 500 12	12
OF 500 14	14
OF 500 16	16
OF 500 18	18
OF 500 20	20
OF 500 22	22
OF 500 24	24
OF 500 26	26
OF 500 28	28
OF 500 30	30
OF 500 32	32
OF 500 34	34

OF 500 36	36
OF 500 38	38
OF 500 40	40
OF 500 42	42
OF 500 44	44
OF 500 46	46
OF 500 48	48
OF 500 50	50
OF 500 52	52
OF 500 54	54
OF 500 56	56
OF 500 58	58
OF 500 60	60



3.5 mm CORTICAL SCREW – SELF TAPPING

Type of bone screw : small screws for small bone and fragment

Fully threaded

Spherical head with small hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 510 08 T	8
OF 510 10 T	10
OF 510 12 T	12
OF 510 14 T	14
OF 510 16 T	16
OF 510 18 T	18
OF 510 20 T	20
OF 510 22 T	22

OF 510 24 T	24
OF 510 26 T	26
OF 510 28 T	28
OF 510 30 T	30
OF 510 32 T	32
OF 510 34 T	34
OF 510 36 T	36
OF 510 38 T	38
OF 510 40 T	40




4.5 mm CORTICAL SCREW – SELF TAPPING

Type of bone screw : small screws for large bone and fragment

Fully threaded

Spherical head with small hexagonal socket

Self tapping with cutting grooves




CODE	Length (mm)
OF 500 10 T	10
OF 500 12 T	12
OF 500 14 T	14
OF 500 16 T	16
OF 500 18 T	18
OF 500 20 T	20
OF 500 22 T	22
OF 500 24 T	24

OF 500 26 T	26
OF 500 28 T	28
OF 500 30 T	30
OF 500 32 T	32
OF 500 34 T	34
OF 500 36 T	36
OF 500 38 T	38
OF 500 40 T	40

3.5 mm CANCELLOUS SCREW – FULL THREADED

Type of bone screws : small screws for small bone fragments

Spherical head with small hexagonal socket



CODE	Length (mm)
OF 512 12	12
OF 512 14	14
OF 512 16	16
OF 512 18	18
OF 512 20	20
OF 512 22	22
OF 512 24	24
OF 512 26	26


OF 512 28	28
OF 512 30	30
OF 512 32	32
OF 512 34	34
OF 512 35	35
OF 512 40	40
OF 512 45	45
OF 512 50	50



3.5 mm CANCELLOUS SCREW – HALF THREAD (1/2 THREADED) – SHORT THREADED

Type of bone screws : small screws for small bone and fragments

Spherical head with small hexagonal socket




CODE	Length (mm)
OF 5121 20	20
OF 5121 22	22
OF 5121 24	24
OF 5121 26	26
OF 5121 28	28
OF 5121 30	30
OF 5121 35	35
OF 5121 40	40
OF 5121 45	45
OF 5121 50	50

4.0 mm CANCELLOUS SCREW – FULL THREADED

Type of bone screws : small screws for small bone and fragments

Spherical head with small hexagonal socket



CODE	Length (mm)
OF 511 10	10
OF 511 12	12
OF 511 14	14
OF 511 16	16
OF 511 18	18
OF 511 20	20
OF 511 22	22
OF 511 24	24


OF 511 26	26
OF 511 28	28
OF 511 30	30
OF 511 35	35
OF 512 35	35
OF 512 40	40
OF 512 45	45
OF 512 50	50



4.0 mm CANCELLOUS SCREW – HALF THREAD (1/2 THREADED) – SHORT THREADED

Type of bone screws : small screws for small bone and fragments

Spherical head with small hexagonal socket



CODE	Length (mm)
OF 5112 20	20
OF 5112 22	22
OF 5112 24	24
OF 5112 26	26
OF 5112 28	28

OF 5112 30	30
OF 5112 35	35
OF 5112 40	40
OF 5112 45	45
OF 5112 50	50


6.5 mm CANCELLOUS SCREW – FULL THREADED

Type of bone screws : large screws for large bone and fragments

Self Tapping with 3 cutting grooves

Spherical head with large hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 501 25	25
OF 501 30	30
OF 501 35	35
OF 501 40	40
OF 501 45	45
OF 501 50	50
OF 501 55	55
OF 501 60	60

OF 501 65	65
OF 501 70	70
OF 501 75	75
OF 501 80	80
OF 501 85	85
OF 501 90	90
OF 501 95	95
OF 501 100	100





6.5 mm CANCELLOUS SCREW – 32 mm THREADED (1/2 THREADED)

Type of bone screws : large screws for large bone and fragments

Self Tapping with 3 cutting grooves

Spherical head with large hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 502 25	25
OF 502 30	30
OF 502 35	35
OF 502 40	40
OF 502 45	45
OF 502 50	50
OF 502 55	55
OF 502 60	60


OF 502 65	65
OF 502 70	70
OF 502 75	75
OF 502 80	80
OF 502 85	85
OF 502 90	90
OF 502 95	95
OF 502 100	100

6.5 mm CANCELLOUS SCREW – 16 mm THREADED (1/3 THREADED)

Type of bone screws : large screws for large bone and fragments

Spherical head with large hexagonal socket

Self Tapping with cutting grooves



CODE	Length (mm)
OF 503 25	25
OF 503 30	30
OF 503 35	35
OF 503 40	40
OF 503 45	45
OF 503 50	50

OF 503 55	55
OF 503 60	60
OF 503 65	65
OF 503 70	70
OF 503 75	75
OF 503 80	80
OF 503 85	85



OF 503 90	90
OF 503 95	95



OF 503 100	100
------------	-----

3.5 mm MALLEOLAR SCREW

Type of bone screws : small screws for small bone and fragments

Design with self a self cutting tip for use in cancellous bone

Spherical head with small hexagonal socket




CODE	Length (mm)
OF 5131 25	25
OF 5131 30	30
OF 5131 35	35
OF 5131 40	40
OF 5131 45	45
OF 5131 50	50
OF 5131 55	55
OF 5131 60	60
OF 5131 65	65
OF 5131 70	70
OF 5131 75	75

4.5 mm MALLEOLAR SCREW

Type of bone screws : small screws for small bone and fragments

Design with self a self cutting tip for use in cancellous bone

Spherical head with large hexagonal socket





CODE	Length (mm)
OF 513 25	25
OF 513 30	30
OF 513 35	35
OF 513 40	40
OF 513 45	45

OF 513 50	50
OF 513 55	55
OF 513 60	60
OF 513 65	65
OF 513 70	70
OF 513 75	75

DHS / DCS LAG SCREW

Used with DHS / DCS Plate



CODE	Length (mm)
OF 502050	50
OF 502055	55
OF 502060	60
OF 502065	65
OF 502070	70
OF 502075	75
OF 502080	80
OF 502085	85
OF 502090	90
OF 502095	95
OF 5020100	100



7. LOCKING SCREW


3.5 mm CORTICAL SCREWS – LOCKING

Type of bone screw : small screws for small bone and fragment

Fully threaded

Spherical head with small hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 510 08 L	8
OF 510 10 L	10
OF 510 12 L	12
OF 510 14 L	14
OF 510 16 L	16

OF 510 18 L	18
OF 510 20 L	20
OF 510 22 L	22
OF 510 24 L	24
OF 510 26 L	26
OF 510 28 L	28

OF 510 30 L	30
OF 510 32 L	32
OF 510 34 L	34
OF 510 36 L	36
OF 510 38 L	38
OF 510 40 L	40


4.5 mm CORTICAL SCREWS – LOCKING

Type of bone screw : small screws for large bone and fragment

Fully threaded

Spherical head with small hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 500 08 L	8
OF 500 10 L	10
OF 500 12 L	12

OF 500 14 L	14
OF 500 16 L	16
OF 500 18 L	18
OF 500 20 L	20



OF 500 22 L	22
OF 500 24 L	24
OF 500 26 L	26
OF 500 28 L	28
OF 500 30 L	30
OF 500 32 L	32
OF 500 34 L	34
OF 500 36 L	36
OF 500 38 L	38
OF 500 40 L	40

OF 500 42 L	42
OF 500 44 L	44
OF 500 46 L	46
OF 500 48 L	48
OF 500 50 L	50
OF 500 52 L	52
OF 500 54 L	54
OF 500 56 L	56
OF 500 58 L	58
OF 500 60 L	60


6.5 mm CANCELLOUS SCREWS – FULL THREAD – LOCKING

Type of bone screws : large screws

Self Tapping with 3 cutting grooves

Spherical head with large hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 501 25 L	25
OF 501 30 L	30
OF 501 35 L	35
OF 501 40 L	40
OF 501 45 L	45
OF 501 50 L	50
OF 501 55 L	55
OF 501 60 L	60

OF 501 65 L	65
OF 501 70 L	70
OF 501 75 L	75
OF 501 80 L	80
OF 501 85 L	85
OF 501 90 L	90
OF 501 95 L	95
OF 501 100 L	100


6.5 mm CANCELLOUS SCREWS – 32 mm THREAD (1/2 THREAD) – LOCKING

Type of bone screws : large screws

Self Tapping with 3 cutting grooves

Spherical head with large hexagonal socket

Self tapping with cutting grooves



CODE	Length (mm)
OF 502 25 L	25
OF 502 30 L	30
OF 502 35 L	35
OF 502 40 L	40
OF 502 45 L	45
OF 502 50 L	50
OF 502 55 L	55
OF 502 60 L	60


OF 502 65 L	65
OF 502 70 L	70
OF 502 75 L	75
OF 502 80 L	80
OF 502 85 L	85
OF 502 90 L	90
OF 502 95 L	95
OF 502 100 L	100

6.5 mm CANCELLOUS SCREWS – 16 mm THREAD (1/3 THREAD) – LOCKING

Type of bone screws : large screws

Spherical head with large hexagonal socket

Self Tapping with cutting grooves



CODE	Length (mm)
OF 503 25 L	25
OF 503 30 L	30
OF 503 35 L	35
OF 503 40 L	40
OF 503 45 L	45
OF 503 50 L	50
OF 503 55 L	55
OF 503 60 L	60

OF 503 65 L	65
OF 503 70 L	70
OF 503 75 L	75
OF 503 80 L	80
OF 503 85 L	85
OF 503 90 L	90
OF 503 95 L	95
OF 503 100 L	100



4 mm CANCELLOUS SCREWS – FULL THREAD – LOCKING

Type of bone screws : small screws for small bone and fragments

Sphericalhead with small hexagonal socket




CODE	Length (mm)
OF 511 10	10
OF 511 12	12
OF 511 14	14
OF 511 16	16
OF 511 18	18
OF 511 20	20
OF 511 22	22
OF 511 24	24
OF 511 26	26
OF 511 28	28
OF 511 30	30
OF 511 35	35

4 mm CANCELLOUS SCREWS – HALF THREAD (1/2 THREAD) – LOCKING

Type of bone screws : small screws for small bone and fragments

Sphericalhead with small hexagonal socket



CODE	Length (mm)
OF 5112 20	20
OF 5112 22	22
OF 5112 24	24
OF 5112 26	26
OF 5112 28	28

OF 5112 30	30
OF 5112 35	35
OF 5112 40	40
OF 5112 45	45
OF 5112 50	50



3.5 mm CANCELLOUS SCREWS – LOCKING

Type of bone screws : small screws for small bone fragments

Spherical head with small hexagonal socket



CODE	Length (mm)
OF 512 12	12
OF 512 14	14
OF 512 16	16
OF 512 18	18
OF 512 20	20
OF 512 22	22
OF 512 24	24
OF 512 26	26
OF 512 28	28
OF 512 30	30




8. PIN

SCHANZ SCREWS

Used for external fixation

With cortical bone thread

Self tapping with cutting grooves



CODE	Length (mm)	Thread Length (L1) (mm)	D Ø (mm)
OF 300 7.5	75	30	3
OF 300 8	80	30	3
OF 300 8.5	85	30	3
OF 300 9	90	30	3
OF 300 10	100	30	3
OF 300 11	110	30	5
OF 300 13	130	30	5
OF 300 15	150	50	5
OF 300 18	180	50	6
OF 300 20	200	50	6
OF 300 22	220	50	6
OF 300 25	250	50	6

STEINMANN PIN

Used for external fixation

With trocar point, slightly rounded up





CODE	Length (mm)	Thread Length (L1) (mm)	D Ø (mm)
OF 310 7.5	75	30	3
OF 310 8	80	30	3
OF 310 8.5	85	30	3
OF 310 9	90	30	3
OF 310 10	100	30	3
OF 310 11	110	30	5
OF 310 13	130	30	5
OF 310 15	150	50	5
OF 310 18	180	50	6
OF 310 20	200	50	6
OF 310 22	220	50	6
OF 310 25	250	50	6



9. ARTHROPLASTY (PROSTHESIS)

SHOULDER STEMMED HEMIARTHROPLASTY – CEMENTED PROSTHESIS


The typical shoulder replacement involves replacing the arthritic joint surfaces with a highly metal ball attached to a stem and should be implanted with bone cement..

Indication :

- For hemiarthroplasty prosthesis of shoulder joint replacement.
- Arthritis that only involves the head of the humerus with a glenoid that has healthy and intact cartilage surface.
- Shoulders with severely weakened bone in the glenoid.
- Some shoulders with severely torn rotator cuff tendons and arthritis

Profile :

Available in various sizes of Humerus head diameter



CODE	Stemmed Length (mm)	Diameter (mm)
OF 4050 35	120	35
OF 4050 36	120	36
OF 4050 37	120	37
OF 4050 38	120	38
OF 4050 39	120	39
OF 4050 40	120	40
OF 4050 41	120	41
OF 4050 42	120	42



ELBOW UNCONSTRAINED PROSTHESIS

2 part prosthesis

Metal to high density with snap fit link

Design consists of – stem for humerus and ulna – resurfacing devices

Fixation with bone cement (bounded to the host bone by bone cement)


Require normal intact ligaments, anterior capsule and appropriate static alignment

Indication :

For total joint replacement prosthesis of elbow reconstruction

Profile :

Part unlinked in a attempt to anatomically duplicates the articular surface of the elbow



CODE	Length (mm)
OF 7000	160



ELBOW SEMI CONSTRAINED PROSTHESIS

Consists of 3 parts prosthesis

Metal to high density articulation connected with locking pin

Fixation with bone cement


Indication :

For total joint replacement prosthesis of elbow reconstruction

Profile :

Available in Right and Left version

Spherical head pin with small hexagonal socket



CODE	Length (mm)	Definition
OF 7100 R	220	Right
OF 7100 L	220	Left



10. EXTERNAL FIXATION

Indications :

1. Fractures and soft tissue care :
 - Open fractures
 - Peri – articular fractures
 - Pediatric fractures
 - Definite treatment
2. Multi – trauma
 - Damage control :
 - i. Pelvic ring injury
 - ii. Long bone fracture temporization
3. Deformity correction :
 - Malunion / nonunion
 - Arthrodesis
 - Limb length
4. Infection
 - Osteomyelitis

1. ROD SYSTEM EXTERNAL FIXATOR

Indication :

Used for tibia and femoral fixation





2. TUBULAR SYSTEM EXTERNAL FIXATOR

Indication :

Used for forearm, metacarpal, and tibia depending on the length of schanz screw and rod used.



3. MINI RAIL EXTERNAL FIXATORY

Stable solution for fractures and for lengthening of small bones. The Mini Rail Fixator is indicated for foot and hand lengthening in adults and for pediatric forearm lengthening.

Mini Rail Single Distractor (Standard)

Indication :

Used for metacarpal, metatarsal, pediatric forearm, phalangeal lengthening and foot surgery.





Mini Rail Articulated X (Vertical)

Indication :

Soft tissue correction in patients with neurological disorders (eg. Cerebral Palsy), arthrodesis of the first carpo – metacarpal joint following soft tissue correction.

This model is mostly used in surgery of the foot (eg. Hallux Valgus)



Mini Rail Articulated Y (Horizontal)

Indication :

Joint stiffness, comminuted intra – articular fractures.





Rotary Single Distractor

Indication :

Used for the treatment of complex distal radius fractures or small bone fracture with single bar distraction.



Rotary Double Distractor

Indication :

Used for forearm, hand and wrist fractures or small bone fractures with double bar distraction.





4. WAGNER

Indication :

Used for tibia and femoral fracture which two connectors in line single bar distractor



5. ILIZAROV APARATUS

Indication :

Used in orthopaedic surgery to lengthen or reshape limb bones, as a limb sparing technique to treat complex and or open bone fractures, and in cases of infected nonunions of bones that are not amenable with other techniques.





11. MISCELLANEOUS IMPLANTS

STAPLES

Indication :

Regular staples which are especially useful in fixing soft tissue to bone, including anterior cruciate reconstruction, also are useful for fragment fixation, tendon and ligament repair and in controlling bone growth.



Osteotomy staples which can be used in high tibial osteotomies to correct varus or valgus deformities, helping to retain cut surfaces in the desired degree of apposition until healing occurs.



Table staples which offer a variety of application including knee ligament repair, arthrodesis of the wrist and ankle, epiphyseal growth arrest, ankle fracture, transplantation of the tibial tubercle, avulsion of the achilles tendon, triple arthrodesis, and proximal humerus fractures.





ENDOBUTTON

Indication :

Used for delivering strength and stiffnes.



WASHER

Indication :

Used to distribute stress under a screw head to prevent thin cortical bone from splitting.

Serrated washers are used to affix avulsed ligaments, small avulsion fractures or comminuted fractures to the remainder of the bone





MANUFACTURING



PT RISA IMPLANTAMA

Jl. Medokan Sawah Timur No. 41 Surabaya

Telp. +62 31 51906646

E mail : osfix.orthopaedics@yahoo.co.id

risa.implantama@gmail.com

Website : www.risaimplantama.com